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6) Anti-corrosion technique: When two different metals, with different chemical properties, come in contact through an electrolyte an electron exchange is established. This process is known as galvanic corrosion which starts oxidation of the weakest metal. In an aluminium autoclave this natural effect can eventually lead to severe damage to the vessel.

To prevent this corrosion, all aluminium autoclave are equipped with a special device called

sacrificial body (Cat. No. PBIB58013). This device is a small disc fixed with a screw to the vessel

body and is made of a metal that is more easily oxidised than aluminium. This means that the

galvanic current will flow through the sacrificial body which will be corroded instead of the autoclave vessel. **If the sacrificial body has been completely worn out, the autoclave vessel will**

begin to corrode, so it is mandatory to periodically check the sacrificial body status and replace it

before it wears out. It is advisable to order some spare sacrificial bodies in order to always run the autoclave with a sacrificial body in it.

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