
## Solar-Grundlagenset 8461200

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09/06 SP

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
### Instruction sheet

09/06 SP

- 1 Holder for metering panel
- 2 Solar metering panel
- 3 Insulating case
- 4 Acrylic plate
- 5 Clamp
- 6 Thermometer

### 1. Safety instructions

The thermometers are sensitive instruments, made from glass. Caution, they are fragile!

- Do not expose to mechanical stress.
- Do not drop.

### 2. Description

The basic solar energy set is an equipment set for experiments on utilisation of solar energy.

The equipment set consists of four solar metering panels of varied surface coating, four thermometers, two holders, two thermal insulating cases and a cover for demonstration of the greenhouse effect.

Four series of measurements, which can be done within approximately 25 min., impressively demonstrate the change in temperature over time and the maximum temperature that can be reached inside the metering panels.

The equipment is delivered in a sturdy cardboard storage box.

### 3. Scope of delivery and technical data

#### 4 Solar metering panels

Material:	Copper
Dimensions:	60 mm x 60 mm
Colour:	1 white, 3 black
Weight:	50 g approx.

#### 2 Insulating cases

Material:	Foam
Dimensions: 1	2 x 12 x 5 cm

#### 1 Acrylic plate

10 x 10 cm

#### 4 Thermometers

-10°C – +100°C

#### 2 Holders for metering panels

#### 2 Clamps

#### 1 Storage box

#### 4. Experimental set-up

##### Recommended accessories:

Halogen lamp with stem, 500 W	8476712
Support base	8614110

- Position the metering panels inside the holders and the insulating cases. The smooth, coloured side should point towards the light source.
- Put the thermometer through the drill hole in the holder or the insulating case, respectively, and into the bottom hole at the rear of the metering panel.
- Mount the acrylic plate with the clamps to one of the insulating cases.
- Position the metering panels at the same distance and angle to the light source.

If the experiment is conducted in direct sunlight, the metering panel can be slanted in the holder to match the angle of the sunbeams.

If the experiment cannot be conducted in direct sunlight, a 500 W halogen lamp has to be used as a substitute.

#### 5. Experimental protocol

- Position the halogen lamp at a distance of approximately 65 cm from the metering panels. The illumination is then as bright as summer sunshine.

- Read thermometer before the experiment, and take a note of the reading.
- Switch on lamp.
- Read the temperature every minute, and enter it into a table.

The various metering panels reach their maximum temperature within different periods of time.

The covered metering panel reaches its maximum temperature after approximately 25 min.

The experiment can be halted once the maximum temperature has been reached.

It can easily be seen (Fig. 1) that black surfaces lead to a much higher rate of temperature increase than white ones.

The thermal insulation of the insulating case prevents energy losses from the rear of the metering panel. Adding the acrylic plate improves utilisation of the radiant power, since the "greenhouse effect" prevents cooling of the meter panel by atmospheric convection and by long-wave heat dissipation from the front. This "greenhouse effect" even compensates for the losses due to absorption by the acrylic plate, which are reflected in the slightly shallower initial slope of the measured curve. The black metering panel, furnished with thermal insulation and a cover, has all the physical attributes of a solar collector panel for a hot water heater.



Fig. 1 Temperature increase in the solar meter panels

1. white, 2. black, 3. black, insulated, 4. black, insulated, with acrylic plate

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