

Instruction manual

VWR® UV-1600PC Spectrophotometer

EU cat. no 634-6001

CE

**UK
CA**

Legal address of Manufacturer

Europe

VWR International bv,
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+3216385011
be.vwr.com

UK Importer

VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
uk.vwr.com

Made in China

Table of Contents

| | |
|------------------------------------------|----|
| Part 1: Spectrophotometer..... | 1 |
| Safety Information..... | 1 |
| Package Contents..... | 2 |
| Unpacking..... | 2 |
| Installation..... | 2 |
| Symbols and Conventions..... | 3 |
| Specifications..... | 3 |
| Overview..... | 4 |
| Description of Buttons and Switches..... | 5 |
| Operational Keys..... | 6 |
| Getting Started..... | 8 |
| Important Guidelines..... | 8 |
| General Operating Instructions..... | 9 |
| Operation..... | 10 |
| Troubleshooting..... | 25 |
| Repair and Maintenance..... | 26 |
| Accessories and Spare Parts..... | 33 |
| Technical Service..... | 34 |
| Warranty..... | 34 |
| Equipment Disposal..... | 35 |
| | |
| Part 2: Software..... | 37 |
| Functions..... | 37 |
| Installation..... | 38 |
| Introduction..... | 44 |
| General Operating Instructions..... | 45 |
| Operation..... | 48 |
| Other Functions..... | 58 |
| Appendix 1..... | 60 |
| Appendix 2..... | 60 |

Part 1: Spectrophotometer

Safety Information

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

VWR recommends against the use of UV-1600PC Spectrophotometer.



- Do not open the device.
- Disconnect the device from the mains supply before carrying out maintenance work or changing the fuses.
- The inside of the device is a high-voltage area Danger!
- Do not use the device if it is damaged, especially if the main power cable way is in any damaged or defective.
- Repairs may only be carried out by the service technicians from your local VWR office and authorized contractual partners.
- The device must be connected to a power outlet that has a protective ground connection.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.



- Do not allow any liquid to enter into the device.
- Do not operate the device in a hazardous location or potentially explosive environment.

Package Contents

| Description | Quantity |
|--------------------------|----------|
| Spectrophotometer | 1PC |
| Glass Cuvette | 4PCS |
| Quartz Cuvette | 2PCS |
| Power Cord (Euro Plug) | 1PC |
| Power Cord (UK Plug) | 1PC |
| Power Cord (CH Plug) | 1PC |
| USB Cable | 1PC |
| CD-ROM | 1PC |
| Quick Manual | 1PC |
| Instruction Manual | 1PC |
| Dust Cover | 1PC |

Unpacking

Open the package, according to carefully check the packaging packing list items, if found inside the packaging are missing or damaged items please contact your local VWR office and authorized contractual partners.

Installation

Placement

Place the instrument on the stable table carefully.

Install printer (Optional)




Check to confirm instrument power switch is turned off, connect the printer's data cable to the Instrument's parallel port.

Link the power cord

Check to confirm instrument power switch is turned off, the power cord plug into two separate power interface and power supply socket apparatus.

Symbols and Conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

| | |
|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
|  | CAUTION This symbol indicates a potential risk and alerts you to proceed with caution |
|  | CAUTION This symbol indicates the presence of high voltage and warns the user to proceed with caution |
|  | CAUTION This symbol indicates risks associated with hot surfaces |

Specifications

| | |
|-----------------------------|--------------------------------------|
| Optical System | Single Beam |
| Wavelength Range | 190-1100nm |
| Band Width | 4nm |
| Stray Light | $\leq 0.05\%T$ @ 220nm & 360nm |
| Photometric Range | 0 to 200%T, -0.3 to 3.0A, 0 to 9999C |
| Wavelength Accuracy | $\pm 0.5nm$ |
| Photometric Accuracy | $\pm 0.5\%T$ or $\pm 0.005A$ @ 1A |
| Stability | 0.002A/h @ 500nm |
| Memory | 200 Results & 200 Standard Curves |

| | |
|-----------------------------|------------------------------------------|
| Language | English, French, German, Spanish |
| Display | 128×64 Dots Matrix LCD |
| Interface | USB, Parallel |
| Measuring Procedure: | Photometry, Quantitation, Kinetics |
| Power Supply | AC 110/220V, 50/60Hz |
| Dimension | 490×360×240 |
| Weight | 14kg |
| Work Environment | 15 to 35°C, 15 to 70% relative humidity |
| Store Environment | -10 to 50°C, 15 to 70% relative humidity |


This instrument is compliant to the European Directives on Low Voltage Directive 2014/35/EU
Electromagnetic compatibility 2014/30/EU
Restriction on use of Hazardous Substances RoHS 2011/65/EU and their amendments.

Overview


UV-1600PC Spectrophotometer used in Chemistry, Pharmaceuticals, Biochemical, Metallurgy, Light Industry, Textile, Material, Environments, Medical, Education and some other fields for Quality Control laboratories.

Description of Buttons and Switches


Front View



Rear View



Operational Keys



Key

Functions

SET

SET Key: Set Parameters

GOTO λ

GOTO λ Key: Set Wavelength



ZERO Key: Blank



PRINT Key: Print measuring result



Numeric Keys: Enter numbers



CLEAR Key: Delete the input value or stored data




RETURN Key: Return to previous interface



START/STOP Key: Start/Stop testing



ENTER Key: Confirm operation



UP, DOWN Keys: Scroll menu/data and set Y scale

Getting Started

The following chart describes the basic operation of the instrument.

Turn On and Self-check

Switch on the power. Then the instrument begins to self-check and 20 minutes' warm up.

Self-check includes the following steps:

√ **Self-testing**

Filter

W Lamp

D2 Lamp

Turn on lamps → Check Sensor → Initialize AD → System position → Get Dark Current → Warm up.

Warm up 20 Minutes,
“Enter” to skip




Important Guidelines

- Reagents and dilution buffers can cause cauterization and other damage to health.
- Samples (nucleic acids, proteins, bacteria cultures) can be infectious and cause serious damage to health.
- During sample preparation, measuring procedures and maintenance and cleaning work, observe all local laboratory safety precautions (e.g. wear protective clothing and gloves, use of disinfectant) regarding the handling of sample material.
- Dispose of measuring solutions and cleaning and disinfectant materials in accordance



with the relevant local laboratory regulations.

General Operating Instructions






Select Application

Main menu, press numeric keys or use the key ,  to choose corresponding menu, then press  to enter into.

Set Wavelength

Press  to set wavelength, use numeric keys to input the values, press  to confirm and go to the point you set, then do blank automatically.

Set Parameters

In different application, press  to set parameters, press ,  to choose or input the values by numeric keys, press  to enter into, press  to return.


Set Auto-cell Holder (optional)

Press the numeric key (1-8) to make corresponding cell position at the light path.

Delete the Input Value

Press  to delete all the characters.


Delete the test result and stored data

Press  to delete the test result or stored data.


Blank

Put the Reference in the light path, press  to do blank.




Measure Samples

Put the samples in the light path, press  to measure.






Print the test results

Press  to print the test results.

Store the Standard Curve

After got the Standard Curve by Marked, press ,  to input the file name and press  to save.

Load the Standard Curve

"Quantitation" interface, Press ,  to select "Load Curve", press ,  to choose the curve you want, press  to open.

Operation

Self-check

Remove all the blocks in the light path and close the lid of the compartment; Switch on the power supply to begin the self-test.

Warm up




After self-test, the instrument goes into pre-warm state. For accurate test, at least 30 minutes of warm up is required.


Check the cuvettes

The cuvettes must be clear and there's no remains of the samples on the surface of it. Only Silicon (Quartz) cuvettes are permitted to be used in the range of UV area.






Photometry

1. Enter into Photometry

Main menu, press numeric key **1** or ,  to choose "Photometry", then press  to go into.



2. Set Photometric Mode

Press  to set photometric mode. Press ,  to choose "Abs.", "T%" or "Energy" and press  to confirm. If user choose "Abs." or "T%". Press  to return.

3. Set Wavelength

Press **GOTO A** to set wavelength, input the value by the numeric keypad followed with **ENTER** pressed to confirm.

4. Blank

Put the Reference in the light path and press **ZERO** to do blank.

5. Measurement samples

Put the sample to be measured in the light path, then the result displays on the screen automatically, press **START STOP** to record.

6. Print Test Results

Press the key **PRINT** to print the Test Results.


```
Photometry
Date&Time: mm-dd-yyyy, hh:mm:ss
Model: UV-1600
Serial No.: UEEXXXXXXX
Firmware Version: 1.5.0
VWR International bvba.
No.  WL.(nm)  Abs  %T
1      500.0  0.000  100.0
2      500.0  0.000  100.0
3      500.0  0.000  100.0
End.
```

Quantitation

1. Enter into Quantitation














Main menu, press **2** or **▲**, **▼** to choose "Quantitation" followed with **ENTER**



pressed to confirm.







2. To establish the standard curve, or call the stored standard curve





Method 1: Use Standard Samples


- 1) Enter into calibration standard samples. Press ,  to select "Standard Curve", then press  to confirm.
- 2) Set Wavelength. Press  to go into wavelength setting interface. Input the value of the wavelength and press  to confirm.
- 3) Blank. Put the Reference Sample in the light path and press  to do blank.
- 4) Set Unit. Press  to set measurement parameters. Press ,  to set "Unit", press  to select, press ,  to choose and press  to select Unit.



5) Setup number of Standard Samples. Press ,  to select "Number", setup number of Standard Samples, press  to confirm.





6) Calibrate Standard Samples. Put the corresponding standard samples in the light path as the screen indicates, input the concentration of the corresponding and press  to measure. Press ,  and  to input curve name and confirm.







Method 2: Input Regression Equation







1) Enter into Coefficient. Press , to select "Coefficient", press to select equation type, then press to confirm.


2) Set test coefficient K. Press to set measurement parameters. Press , to select "Coefficient K", input test coefficient K, press to confirm.





3) Set test coefficient B. Press ,  to select "Coefficient B", input test coefficient K, press  to confirm.









4) Set Unit. Press ,  to set "Unit", press  to select, press ,  to choose and press  to select Unit.



5) Set Wavelength. Press  to go into wavelength setting interface. Input the value of the wavelength and press  to confirm.


Load the Stored Curve

Press ,  to choose "Select Curve", press ,  to choose the cursor, press  to confirm.



3. Return the sample measurement interface


Press  to return the Quantitative Test interface.




4. Blank

Put the reference in the light path, press  to do blank.

5. Measure Samples

Place the sample to be tested in the light path, press  to measure. Then the test result will display in the data sheet. Repeat this step to finish measuring all the samples.

6. Print Test Results


Press the key  to print the Test Results.

```
Quantitation
Date&Time: mm-dd-yyyy, hh:mm:ss
Model: UV-1600PC
Serial No.: UEEXXXXXXX
Firmware Version: 1.5.0
VWR International bvba.
C=1.000*A+1.000
r=1.0000
No.  WL.(nm)  Abs  C(mg/L)
1     500.0  0.120  1.120
2     500.0  0.127  1.127
3     500.0  0.121  1.121
End.
```

Kinetics

1. Enter into Kinetics

Main menu, press  or ,  to select "Kinetics" and press  to confirm.






2. Set Wavelength

Press  to set wavelength, input the value of the wavelength by numeric keypad and

press **ENTER** to confirm.

3. Setup Parameters

Press **SET** to set measurement parameters. Press ,  to select values of "Time Intervals" , "Total Time", "Test Mode", "Upper Limited" and "Lower Limited". Press **ENTER** to confirm.




4. Blank


Put the Reference Sample in the light path, press **ZERO** to do blank.

5. Measure Samples

Put the sample in the light path and press **START STOP** to enter the measurement interface, press **START STOP** to begin the test, repress it to stop, press **RETURN** to cancel.





6. Print Test Results

Press the key  to print the Test Results.


```
Kinetics
Date&Time: mm-dd-yyyy, hh:mm:ss
Model: UV-1600PC
Serial No.: UEEXXXXXXX
Firmware Version: 1.5.0
VWR International bvba.
WL=500nm
Total Time: 5s
Interval: 1s
No.  Time(s)  Abs    %T
1      1.0  0.000  100.0
2      2.0  0.000  100.0
3      3.0  0.000  100.0
4      4.0  0.000  100.0
5      5.0  0.000  100.0
End.
```


7. View Values

After scan finished, press  and press ,  to select "Curve Data", press  to list all the values..








Utility

Main menu, press **4** or use ,  to select "System Utility" and press **ENTER** to confirm.








Turn On/Off D2 Lamp

Press ,  to choose "D2 Lamp", then press **ENTER** to enter into. Press ,  to choose "On" or "Off", press **ENTER** to turn on/off.








Turn On/Off W Lamp

Press ,  to choose "W Lamp", then press **ENTER** to enter into. Press ,  to choose "On" or "Off", .press **ENTER** to turn on/off.






Set Date & Time


Press ,  to choose "Date & Time", then press  to enter into. Input hh:mm:ss and yy-mm-dd, press  to confirm.






Get Dark Current


Keep the light path without anything blocking, press ,  to choose "Dark Current", then press  to resample Dark Current.

Note: During the course, open the lid of the compartment is prohibited.







Reset Wavelength

Keep the light path without anything blocking, press ,  to choose "Reset WL", then press  to reset wavelength.







Lamp Life

Press ,  to choose "Lamp Life", then press  to enter into. Select "Display Lamp Life" to view the W Lamp and D2 Lamp has been used time. Select "Reset D2 Life" or "Reset W Life" to reset the time.








Load Default Parameter

Press ,  to choose "Load Default", then press  to enter into. Select "Yes" to load the parameters to factory setting and the Instrument will restart.




Lamp Change

Press ,  to choose "Lamp Change", then press  to enter into. Input switch wavelength(325—375nm)and press  to confirm.



About Version

Press ,  to choose "Version", press  to view version information, press any key to return.



Troubleshooting

Review the information in the table below to troubleshoot operating problems.

| Problem | Cause | Solution |
|-------------------------|-----------------------------------------|--------------------------|
| Power on, no response | Power cord connection is not reliable | Improve connectivity |
| | Fuse burning | Replace fuse |
| Measurement uncertainty | Warm up is not enough | Warm up more time |
| | Glass cuvettes used in UV Range | Use quartz cuvettes |
| | Sample is not Stable | Improve the sample |
| | The concentration of sample is too high | Diluted sample |
| | Power Supply Voltage Low or not Stable | Improve the Power Supply |
| | Lamp damage or lamp life maturity | Replace lamp |
| Dark Current Error | The lid of the compartment is | Close the lid, restart |

| | | |
|---------------------------------------------------------------------------------------|--------------------------------|--------------------------------------|
| when self-check | open during self-check | |
| System Calibrate Failed | Something block the Light path | Remove it, calibrate again |
| Power on, back light is OK, but nothing display on the screen or display is not clear | Display Contrast problem | Adjust the contrast potentiometer |
| Measurements inaccurate | Cuvettes were contaminated | Clean cuvettes |
| | Samples were contaminated | Improve samples |
| | Worse matching of the cuvettes | Improve the matching of the cuvettes |
| | Dark current error | Resample dark current |

Repair and Maintenance

Daily Maintain

Check the compartment

After measurement, the cuvettes with sample solutions should be taken out of the compartment in time. Or the volatilization of the solution would make the mirror go moldy. Users must pay more attention to the corrosive sample and liquid easy to volatilize. Any solution remains in the compartment should be wipe off immediately.

Surface clean

The cover of the instrument is with paint. Please use wet towel to wipe off the drips on the surface immediately. Organic solution is forbidden to be used to clean the cover. Please wipe off the dirt on the cover timely.

Clean the cuvettes

After every test or after a solution change, the cuvettes should be cleaned carefully, or the remains on the surface would cause measuring error.

Check Lamp

Check the W lamp

In "Photometry" mode, set the test parameters as follow:

WL.: 500nm

Test Mode: Energy

Gain: 1

Check the energy, you need replace W lamp as following two cases:

Energy < 20 W lamp is damaged

Energy < 5000 W lamp energy is too low

Check the D2 lamp

In "Photometry" mode, set the test parameters as follow:

WL.: 200nm

Test Mode: Energy

Gain: 1

Check the energy, you need replace D2 lamp as following two cases:

Energy < 20 W lamp is damaged

Energy < 1000 W lamp energy is too low

Spare Parts Replacement

Replace the fuse



Danger! Be sure to switch off the power and unplug the socket before replacement!

1. Tools preparation


Prepare a 3x75 Flat Blade screwdriver.

2. Switch Off the power supply

Switch off the power supply, and unplug the socket.


3. Take out the Fuse Seat

Take out the fuse seat by the screwdriver.



4. Replace a new fuse

Pick out the spare fuse (3.15A/250V) and replace it to the working position.



5. Reset the fuse seat

Replace the fuse seat in the power socket.

6. Switch on the power

Plug the socket and switch on the power.

Replace lamps



Hot! Wait 20 minutes before open the lamp chamber after power off to avoid scald!

1. Tools preparation


Prepare a 6×150mm Cross Blade screwdriver and a pair of glove.

2. Power Off

Switch off the power supply and unplug the socket.


3. Open the cover

Unscrew the 4 screws indicated(Each side with 2 screws) and remove the cover.




4. Open the cover of the light chamber

Unscrew the 2 screws on the light chamber cover and remove it.



5. Replace the D2 lamp



Unscrew the 2 screws on the D2 Flange (No.1), unplug the connector in the power board(No. 2) and remove the D2 lamp. Draw on the cotton glove and replace a new lamp. Fix the 2 screws and plug the connector again.



6. Replace W lamp

The Tungsten lamp is equipped with a blue-grey silicon coating by manufacturer. This coating is only a transport safety device. It can be removed with the first exchange of lamp.



Pull out the defected W lamp and draw on the cotton glove. Insert the new W lamp as deep as possible on the lamp seat. Be sure to keep the filament in the same direction as the old one face.



Filament

7. Adjust the position of the W lamp

Switch on the power(the Switch Mirror should be placed to the position as indicates). Observe the entrance facular, and it should in the center of the entrance hole. If the facular deviate to Left or Right, then loosen the No.1 screws in Fig. 5-8 and move the lamp seat to Left or Right until it focus on the center of the slot. Then fix the screws. If the facular deviate to Up and Down, then loosen the No.2 screws and move the lamp seat Up and Down until the facular focus on the center of the slot. Then fix the No. 2 screws again.



8. Finish

Reset the cover of the light chamber and fix the screws. Reset the cover of the instrument and fix the screws. Recover the Pole in the compartment, then the course finished.

Replace the Battery



Be sure to switch off the power supply and unplug the socket before open the Bottom Cover!

1. Prepare the tools

Prepare a 6×150mm Cross Blade Screwdriver.

2. Switch off the power supply

Switch off the power supply and unplug the socket


3. Open the Bottom cover plate

Unscrew all the screws indicated then remove the bottom plate.



4. Replace the Battery

Pick out the old battery and replace a new one.



5. Finish

Recover the bottom plate and fix all the screws, then the course finishes.

Accessories and Spare Parts

| Description | Quantity | Cat. No. |
|----------------------------------------------|----------|----------|
| CELL HOLDER, 4-CELL, 10MM | 1PC | 634-6003 |
| CELL HOLDER, 4-CELL,10 TO 50MM | 1PC | 634-6004 |
| CELL HOLDER, 4-CELL, 10 TO 100MM | 1PC | 634-6005 |
| CELL HOLDER, FOR CYLINDRICAL CELL | 1PC | 634-6006 |
| CELL HOLDER, WATER-JACKETED, 1 CELL, 10MM | 1PC | 634-6007 |
| CELL HOLDER FOR MICRO CELLS,BEAM HEIGHT:15MM | 1PC | 634-0687 |
| CELL HOLDER FOR TEST TUBES | 1PC | 634-6009 |
| CELL HOLDER,8-POSITION AUTO CELL CHANGER | 1PC | 634-6010 |
| CELL HOLDER, SOLID SAMPLE | 1PC | 634-6011 |
| CELL HOLDER, WATER-JACKED, 4 CELL, 10MM | 1PC | 634-6012 |
| CUVETTE, SQUARE.GLASS,10MM | 4PCS | 634-6013 |
| CUVETTE, SQUARE.GLASS,20MM | 2PCS | 634-6014 |
| CUVETTE, SQUARE.GLASS,30MM | 2PCS | 634-6015 |
| CUVETTE, SQUARE.GLASS,50MM | 2PCS | 634-6016 |
| CUVETTE, SQUARE.GLASS,100MM | 1PC | 634-6017 |
| CUVETTE, SQUARE.QUARTZ,10MM | 2PCS | 634-6018 |
| CUVETTE, SQUARE.QUARTZ,20MM | 2PCS | 634-6019 |
| CUVETTE, SQUARE.QUARTZ,30MM | 2PCS | 634-6020 |
| CUVETTE, SQUARE.QUARTZ,50MM | 2PCS | 634-6021 |
| CUVETTE, SQUARE.QUARTZ,100MM | 1PC | 634-6022 |
| CELL,QUARTZ, 100UL, 10MM, BEAM HEIGHT:15MM | 1PC | 634-0688 |
| CELL,QUARTZ, 200UL, 10MM, BEAM HEIGHT:15MM | 1PC | 634-0689 |
| CELL,QUARTZ, 500UL, 10MM | 1PC | 634-6025 |
| FLOW CELL, 5 MM, GLASS, BEAM HEIGHT:15MM | 1PC | 634-0690 |
| FLOW CELL, 10 MM, GLASS, BEAM HEIGHT:15MM | 1PC | 634-0691 |
| FLOW CELL, 20 MM, GLASS, BEAM HEIGHT:15MM | 1PC | 634-0692 |
| FLOW CELL, 30 MM, GLASS, BEAM HEIGHT:15MM | 1PC | 634-0693 |
| FLOW CELL, 5 MM, QUARTZ, BEAM HEIGHT:15MM | 1PC | 634-0694 |

| | | |
|----------------------------------------------------------|-----|----------|
| FLOW CELL, 10 MM, QUARTZ, BEAM HEIGHT:15MM | 1PC | 634-0695 |
| FLOW CELL, 20 MM, QUARTZ, BEAM HEIGHT:15MM | 1PC | 634-0696 |
| FLOW CELL, 30 MM, QUARTZ, BEAM HEIGHT:15MM | 1PC | 634-0697 |
| PELTIER UNIT, BEAM HEIGHT:15MM | 1PC | 634-0698 |
| SIPPER UNIT WITHOUT TEMP. CONTROL, BEAM HEIGHT:15MM | 1PC | 634-0699 |
| SIPPER UNIT WITH PELTIER TEMP. CONTROL, BEAM HEIGHT:15MM | 1PC | 634-0763 |
| LAMP, HALOGEN, 12V20W | 1PC | 634-6037 |
| LAMP, DEUTERIUM | 1PC | 634-6038 |
| PRINTER, THERMAL PRINTER | 1PC | 634-6039 |
| THERMAL PAPER | 1PC | 634-6043 |
| FUSE, 3.15A/250V | 1PC | 634-0651 |
| FUSE, 500MA/250V | 1PC | 634-0652 |
| BATTERY, 3V, CR2032 | 1PC | 634-0653 |

Technical Service

Web Resources

Visit the VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

Contact us For information or technical assistance contact your local VWR representative or visit. www.vwr.com.

Warranty

VWR International warrants that this product will be free from defects in material and

workmanship for a period of two (2) years from date of delivery except the lamps. Lamps have a warranty of 1000 h lamp usage time or 6 months max. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non-performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the Product in its local environment. VWR will not be held liable for any related omission or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment Disposal



This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment. For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you!

Part 2: Software

Functions

We'll introduce the main functions of M.Wave Professional in this chapter.

Quantitative Analysis

- Provide 2 methods to establish a Standard Curve;
- Up to 20 standard samples to establish a new curve or input the coefficients to make a standard curve;
- There are three methods, which are Linear zero , Linear and Square to fit the standard curve.

Kinetics Analysis

- We can choose the Time Intervals(0.5, 1.0, 2.0, 5.0, 10.0, 30.0 or 60.0 Seconds);
- You can choose different Photometric Mode to display the curve (Transmittance-Time & Absorbance-Time).

Spectrum Scan

- You can choose the Scan Intervals(0.1, 0.2, 0.5, 1.0, 2.0 and 5nm);
- You can choose the Photometric Mode to display the spectrum (Wavelength-Transmittance, Wavelength-Absorbance);
- Automatically list spectrum peaks;
- Mathematics and smooth.

Multi-wavelength Analysis

- Up to 20 wavelengths to measure a sample.

DNA/Protein Analysis

- We provide 2 methods;
- You can input the coefficients by yourself.

Energy Scan

- You can switch the light source point(W Lamp, D2 Lamp or automatically

switch);

- You can choose the scanning intervals(0.1,0.2,0.5,1.0,2.0 and 5.0);
- You can set the gain(1 to 8 times);
- Automatically list spectrum peaks;

Installation


We will show you how to install M.Wave Professional on your PC in this chapter.

PC System Requirements


- Pentium or above PC;
- CD-ROM;
- USB Ports.
- 512MB Memory(1GB or Above is strongly recommended);
- 50MB or above hard disc Space;
- Microsoft Windows XP/Vista/7/8.

Install M.Wave Professional


1. Put M.Wave Professional disc in the CD-ROM;
2. Double Click to open the CD-ROM, and then, double click **Setup.exe** which is under the root directory of CD to start installation , "Notice" form appear, click **OK**;




3. "Welcome" form appear, click **Next**;




4. Input user's information, click **Next**;




5. Choose install path, then click **Next**;




6. Choose Shortcut file folder, then click **Next**, all the setting information will be displayed on the screen, click **Next** to copy files to Hard disc;




7. After finished copy all files of M.Wave Professional, it will start install the USB Drive. Click **Next**;




8. Select "I accept this agreement". Click **Next** to copy files to PC;




9. Click **Finish** to finish install USB Drive;



10. "Notice" form appear, click **OK**;




11. Click **Finish** to finish all the setup.



Uninstall M.Wave Professional

1. Start→All Files→M.Wave Professional→Uninstall M.Wave Professional 2.0, click **Next** to uninstall.
2. After all the files were removed, click **Finish** to end.




Run M.Wave Professional




- Double click the icon on desktop.
- Start→All Files→M.Wave Professional→M.Wave Professional 2.0.

Set up Communication Port

Be sure that the USB line which is connecting PC with instrument is working and the spectrophotometer is at the standby state interface.




After M.Wave Professional being started, then click **Main Menu→Options**, a dialog box of **Options** will display in the screen, click **Search** to automatically search the instrument's communication port, after get the communication port, click **OK** to save setting. If you choose the option of Connect to Spectrophotometer When Start, the instrument will connect to the PC automatically when you run M.Wave Professional next time.




Set User's Information

Click **Main Menu**→**Options**, pop-up **Options** window, click **User's Information** tag, choose it and input relative user information, click **OK** to save settings, all this information will appear in your testing report.




Set Data Format

Click **Main Menu**→**Options**, a dialog box of **Options** will display in the screen, click **Data Format** tag, choose the displaying format, and click **OK** to save setting.




Set Interface Style

Click **Main Menu**→**Options**, a dialog box of **Options** will display in the screen, and click **Interface** tag. You can custom the style and color schemes of the interface as you like. Click **OK** to save setting.



Connect to PC


After enter the Main Interface, click  in the Tool Bar to connect with PC. After connection, the icon is on the pitched on state, you can release PC if you click the icon again.

Introduction

We will introduce the M.Wave Professional in this chapter.

Main Interface



This is the main interface after start.




Menus Bar and Tools Bar

Menus Bar and Tools Bar provide three different ways to operate this software.


- You can finish all the operation through Keyboard or Mouse.
- Most of the functions can be finished through the Shortcut Bar.
- Most of the functions can be finished through the menu by Right-hand button.







General Operating Instructions

Set Wavelength

Click the icon  on the Tools Bar to set wavelength. Input the value of wavelength in the frame of **Set Wavelength**, click **Goto** to confirm. Then the instrument will move the wavelength to the point you want, then it blanks automatically.



Blank

Put reference into the light path, click the icon  in the Tool Bar.


New a Testing

Main Menu→File→New, select a testing to click.

Lists of the file's type

- Quantity Analysis File *.qua
- Kinetics Analysis File *.kin
- Spectrum Scan File *.wls
- Multi-wavelength File *.mti
- DNA/Protein Analysis File *.dna
- Energy Scan File *.ens
- Standard Curve *.std


Test a Sample

Put sample into the light path, click the icon  in the Tool Bar.


Set Test Parameters

Click the button of  to set the parameters.

Modify a Sample

1. Choose the data frame you want to modify or set the curve as current curve;
2. Pull the unknown sample in light path;
3. Click the icon  in Tools Bar to re-measure it, and save the current result to substitute the original one.


Delete a Test Result

1. Choose the data frame you want to modify or set the curve as current curve;
2. Click the icon  in the Tools Bar to delete the data.


Name a Sample

1. Select the line you want to rename.
2. Double-click **Sample Name**, input the sample's name, and press **Enter** on the keyboard to confirm.

Save a Testing

Click the icon  in Tools Bar, input the file name and click **Save** to confirm in the dialog box.

Open a Testing

Click the icon  in Tools Bar, choose the file name you need in the dialog box, then click **Open** to confirm.

Print Test Report

Click the icon  in Tools Bar, then set Print Parameters in the dialog box, press **Print** to confirm.

Export Data to Excel (the Microsoft Excel should be installed)

Click **Main Menu**→**File**→**Export to Microsoft Excel**, Excel will be started up, and data will be stored in it.



Export Curve as bmp Picture

Click **Main Menu**→**File**→**Export Bitmap**, input the file name and click save to confirm the dialogue box.


Switch Display Mode


Click the icon **A** or **T** in the Tools Bar to switch the display mode (Abs. or Transmit).

Zoom Selected Area

Click the icon  in Tools Bar, choose the area you want to magnify, then you'll see the magnified one, the Mouse's cursor turns into the form of "+". Click the icon  again to cancel the current status.

Modify Scale

Click the icon  on Tools Bar to display coordinate as you like.



Auto Fitting

Click the icon on Tools Bar to change the scale to appropriate value.



Select as Current Curve



Click **Main Menu**→**Curve Operate**→Choose the curve as you like.

Change the Color of Current Spectrum






Click **Main Menu**→**Curve Operate**→Choose the color as you like.





Search Peaks

Click the icon  in Tools Bar to search the peaks, they labeled automatically on the curve, and the peaks will be listed in the Right Frame. You can set the Peak height by clicking  icon.



Mathematics Function

Click the icons     or  in the Tools Bar, a new dialog box will display on the screen respectively. You can add, subtract, multiply or divide two spectrums.



Spectrum Smooth

We use this function to wipe off the interferer of noise during the course of scan, so as to make the spectrum look much smoother. Click **Main Menu**→**Spectrum**→**Spectrum Smooth**.


Operation

We will introduce you how to use M.Wave Professional.


Quantitative Analysis

M.Wave Professional use Standard Curve to test the sample's concentration at a fixed wavelength point.

1. Click the icon  in the Tools Bar shortcut to create a new Quantitative Analysis.



2. Click the icon  in the Tools Bar Shortcut to set the parameters of Quantitative Analysis.



The 'Quantitative Analysis Setup' dialog box is shown. It has several sections: 'Method' with radio buttons for 'Single Wavelength', 'Double Wavelength', and 'Three Wavelength'; 'Unit' set to 'µg/ml'; 'Conc.' set to 'Concentration'; 'Fitting' set to 'Linear with 0'; 'WL.' section with 'WL. 1' set to '300.8', 'WL. 2' set to '550', and 'WL. 3' set to '600'; 'Coefficient' section with 'k0' set to '0', 'k1' set to '1.0000', and 'k2' set to '0'; and 'Standard Sample Quantity' set to '3'. There are also fields for 'Conc. 1' through 'Conc. 20'.

- Choose test method;
- Input test wavelength in **WL** Frame;
- Choose the Concentration Unit ;
- Choose curve fitting method;

Create Standard Curve

Two methods are under your choice to create a new curve.

Method 1: Coefficient

- (1) Click Coefficient option;
- (2) Click Fit Method to choose Fit Curve method;
- (3) Input the curve equation's coefficients in corresponding frames;
- (4) Click **OK** to finish setting.

Coefficient

k0

k1

k2



Method 2: Use Standard Samples to create a new curve

- (1) Click **Standard Sampling** option;
- (2) Click **Standard Sample Quantity** to choose the quantity(≤ 20);
- (3) Input the sample concentration in corresponding Frame;




Conc.


Standard Sample Quantity



| | | | |
|----------|------------------------------------|----------|---------------------------------|
| Conc. 1 | <input type="text" value="1.000"/> | Conc. 11 | <input type="text" value="11"/> |
| Conc. 2 | <input type="text" value="2.000"/> | Conc. 12 | <input type="text" value="12"/> |
| Conc. 3 | <input type="text" value="3.000"/> | Conc. 13 | <input type="text" value="13"/> |
| Conc. 4 | <input type="text" value="4"/> | Conc. 14 | <input type="text" value="14"/> |
| Conc. 5 | <input type="text" value="5"/> | Conc. 15 | <input type="text" value="15"/> |
| Conc. 6 | <input type="text" value="6"/> | Conc. 16 | <input type="text" value="16"/> |
| Conc. 7 | <input type="text" value="7"/> | Conc. 17 | <input type="text" value="17"/> |
| Conc. 8 | <input type="text" value="8"/> | Conc. 18 | <input type="text" value="18"/> |
| Conc. 9 | <input type="text" value="9"/> | Conc. 19 | <input type="text" value="19"/> |
| Conc. 10 | <input type="text" value="10"/> | Conc. 20 | <input type="text" value="20"/> |


- (4) Click **OK** to finish setting;
- (5) Pull the Reference into the light path and click the icon , then the system will move the wavelength to the one you just set, it blanks then;
- (6) Put **1# standard sample** into the light path and move the cursor on the first frame of Abs. The sample Absorbance value will appear in this frame after click  in the shortcut toolbar;
- (7) Measure other standard samples in the same way as step (6) shows. The curve displays automatically when finish;

| ID | Abs. | Conc. (ug/ml) |
|----|--------|---------------|
| 1 | 0.0352 | 1.00 |
| 2 | 0.1887 | 2.00 |
| 3 | 0.3075 | 3.00 |

You can modify the curve's coordinates by Clicking the icon . You can also save the curve by clicking , and you can reload it by clicking  next time.





3. Pull the unknown into the light path (if it is Automatic 8-cell Holder in the instruments, click  to set the location of the relative cell first), Then click the icon, click  in the shortcut toolbar to test, the test result will display in the data sheet.




Kinetics Analysis


We will introduce how to measure a sample's change over a selected period of time.


1. Click the icon  to begin a new test of Kinetics;



2. Click the button of  to set the parameters;





3. Choose Photometric Display Mode, input the Wavelength, Display Time, Total Time, Intervals, Responsible Mode and Coordinates Values;
4. Click OK to finish and exit setting;
5. Pull the reference into the light path and click icon  in the Tool Bar to do blank.
6. Pull the sample into the light path and click the button of to begin the test. You can click icon in the Tools Bar to cancel at any time.





Spectrum Scan


We'll introduce how to create the spectrum while using Spectrum Scan function.


1. Click the icon  on the Tools Bar to start a new Spectrum scan;






2. Click the icon  in the Tool Bar to set the parameters of Wavelength Scan;



3. Choose the Photometric Mode, Wavelength Range for Scanning, Intervals, Repeat Times, Response Mode, and Coordinate Range;
4. Click **OK** to finish and exit setting;
5. Put reference into the light path, click the icon  in the Tool Bar, then the dialog box of **Baseline** displays in the screen, click **Scan** to scan baseline, click **Cancel** to interrupt scanning and exit the course;




6. Put the unknown sample in the light path, click the icon  in the Tools Bar to start scanning, click  to interrupt and **cancel** scanning.



Multi-wavelength Analysis

We'll introduce how to measure a sample in different wavelengths (≤ 20).

- Click the icon  on the Tools Bar, a new form appears on the screen;





- Click the icon  on the Tools Bar to set the parameters;


Multi Wavelength Analysis

Setup

WL No. 3

| | |
|------------------------------------------------------------------|------------------------------------------------------------------|
| WL.1(nm) <input style="width: 80%;" type="text" value="500.0"/> | WL.11(nm) <input style="width: 80%;" type="text" value="600.0"/> |
| WL.2(nm) <input style="width: 80%;" type="text" value="510.0"/> | WL.12(nm) <input style="width: 80%;" type="text" value="610.0"/> |
| WL.3(nm) <input style="width: 80%;" type="text" value="520.0"/> | WL.13(nm) <input style="width: 80%;" type="text" value="620.0"/> |
| WL.4(nm) <input style="width: 80%;" type="text" value="530.0"/> | WL.14(nm) <input style="width: 80%;" type="text" value="630.0"/> |
| WL.5(nm) <input style="width: 80%;" type="text" value="540.0"/> | WL.15(nm) <input style="width: 80%;" type="text" value="640.0"/> |
| WL.6(nm) <input style="width: 80%;" type="text" value="550.0"/> | WL.16(nm) <input style="width: 80%;" type="text" value="650.0"/> |
| WL.7(nm) <input style="width: 80%;" type="text" value="560.0"/> | WL.17(nm) <input style="width: 80%;" type="text" value="660.0"/> |
| WL.8(nm) <input style="width: 80%;" type="text" value="570.0"/> | WL.18(nm) <input style="width: 80%;" type="text" value="670.0"/> |
| WL.9(nm) <input style="width: 80%;" type="text" value="580.0"/> | WL.19(nm) <input style="width: 80%;" type="text" value="680.0"/> |
| WL.10(nm) <input style="width: 80%;" type="text" value="590.0"/> | WL.20(nm) <input style="width: 80%;" type="text" value="690.0"/> |


3. Choose the quantity of the wavelength, you need input the value of each wavelength;
4. Click **OK** to finish and exit setting;
5. Pull the reference into the light path, click the icon  in the Tools Bar to do blank;
6. Pull the sample into the light path, click  to start test. And the test result will be listed in the data sheet.




| ID | Sample Name | Abs. (200.0nm) | Abs. (210.0nm) | Abs. (220.0nm) |
|----|-------------|----------------|----------------|----------------|
| 1 | | -0.0002 | 0.0002 | 0.0009 |
| 2 | | 0.0012 | 0.0009 | 0.0004 |
| 3 | | 0.0070 | 0.0008 | 0.0059 |
| 4 | | 0.0026 | 0.0022 | 0.0015 |
| 5 | | 0.0020 | 0.0013 | 0.0011 |
| 6 | | -0.0020 | -0.0019 | -0.0008 |
| 7 | | -0.0002 | 0.0000 | -0.0001 |
| 8 | | -0.0002 | 0.0001 | 0.0000 |

DNA/Protein Analysis


We will introduce how to operate a DNA/Protein Analysis.



1. Click the icon  on the Tools Bar to begin a DNA/Protein Analysis;



| ID | Sample Name | Abs. (260.0nm) | Abs. (280.0nm) | Abs. (320.0nm) | DNA Conc. (µg/ml) | Protein Conc. (µg/ml) | Ratio |
|----|-------------|----------------|----------------|----------------|-------------------|-----------------------|-------|
| | | -0.0007 | | 100.17% | | | |

2. Click the icon  to set parameters;




3. Choose the measure method first. We provide you 2 methods to measure, and you can also definite the coefficients and wavelength by yourself;
4. Click **OK** to finish and exit the setting;
5. Pull the reference into the light path, click  in the Tools Bar to do blank;
6. Pull the sample in the light path, click the icon  to measure. The result will display in the data sheet.


| ID | Sample Name | Abs. (260.0nm) | Abs. (280.0nm) | DNA Conc. (ug/ml) | Protein Conc. (ug/ml) | Ratio |
|----|-------------|----------------|----------------|-------------------|-----------------------|-------|
| 1 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 2 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 3 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 4 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 5 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 6 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 7 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |
| 8 | | 3.0000 | 3.0000 | 80.700 | 2384.100 | 1.000 |


Energy Scan

We will introduce how to scan the variability of energy in the range of wavelength.

1. Click **Main Menu**→**File**→**New**→**Energy Scan** to build Energy Scanning;



2. Click the icon  on the Tools Bar to set parameter;






3. Select Lamp Switch Mode, Wavelength Range, Intervals, Repeat Times and Coordinate Range;



- **Automatic Switch Mode**---The instrument can choose D2 Lamp or W Lamp according to setting;
- **Deuterium Lamp Mode**---Using Deuterium Lamp in the whole wavelength;
- **Tungsten Lamp Mode**---Using Tungsten Lamp in whole wavelength.

4. Click **OK** to finish and exit setting;

5. Pull the unknown sample in the light path, click the icon  in Tools Bar to start scanning, click  to interrupt scanning and exit the course;



Other Functions

We will introduce the additional function about M.Wave Professional in this chapter.

Turn On/Off W Lamp


Click **Main Menu**→**Device**→**Turn On/Off W Lamp** to open or close W Lamp.

Turn On/Off D2 Lamp

Click **Main Menu**→**Device**→**Turn On/Off D2 Lamp** to open or close D2 Lamp.

Setup Switch Point of Light Source

Click **Main Menu**→**Device**→**Setup Switch Point of Light Source**, input the value of wavelength, and then click **OK** to confirm.




Get Dark Current

Click **Main Menu**→**Device**→**Get Dark Current**, the system will get the dark current and the new data will replace the old one.



Scan System Baseline

Click **Main Menu**→**Device**→**System Baseline**, the dialogue box will display in the screen, you can click **Open** to choose the storage base line to test, be sure that nothing in the light path, you can click **Scan** to build system base line. You can click **Cancel** to interrupt and exit the course. You can choose **Save System Base Line** option; a dialogue box will display in the screen and click **Save** to save system base line for your after use.



You needn't to build system base line normally, only when it is a very long time that the instrument wasn't be used or the environment changed much.

Automatic 8-Position Changer (Accessories)

Click the icon  ...  to pull the correspond cell in the light path.

Cell Calibration (Need install auto 8 cells Changer)

Click **Main Menu**→**Device**→**Cell Calibration**,this operation for matching difference of all the cells selected.

Appendix 1

Methods of Quantitative Analysis

Single Wavelength Method: $Abs.=A_1$

Double Wavelengths Method: $Abs.=m \cdot A_1 - n \cdot A_2$

Three Wavelengths Method: $Abs.=A_1 - (WL_1 - WL_2) \cdot (A_2 - A_3) / (WL_2 - WL_3) - A_3$

Appendix 2

Methods of DNA/Protein Analysis

Method 1: $C_{DNA} = (A_1 - A_{ref}) \cdot f_1 - (A_2 - A_{ref}) \cdot f_2$
 $C_{Protein} = (A_2 - A_{ref}) \cdot f_3 - (A_1 - A_{ref}) \cdot f_4$
 $Ratio = (A_1 - A_{ref}) / (A_2 - A_{ref})$

$A_1 = A_{260nm}$, $A_2 = A_{280nm}$, $A_{ref} = A_{320nm}$ (Optional)
 $f_1 = 62.9$, $f_2 = 36.0$, $f_3 = 1552$, $f_4 = 757.3$

Method 2: $C_{DNA} = (A_1 - A_{ref}) \cdot f_1 - (A_2 - A_{ref}) \cdot f_2$
 $C_{Protein} = (A_2 - A_{ref}) \cdot f_3 - (A_1 - A_{ref}) \cdot f_4$
 $Ratio = (A_1 - A_{ref}) / (A_2 - A_{ref})$

$A_1 = A_{260nm}$, $A_2 = A_{230nm}$, $A_{ref} = A_{320nm}$ (Optional)
 $f_1 = 49.1$, $f_2 = 3.48$, $f_3 = 183$, $f_4 = 75.8$

Local VWR offices in Europe and Asia Pacific

Austria

VWR International GmbH
Graumannsgasse 7
1150 Vienna
Tel.: +43 1 97 002 0
Email: info.at@vwr.com

Belgium

VWR International BV
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 016 385 011
Email: vwr.be@vwr.com

China

VWR (Shanghai) Co., Ltd
Bld.No.1, No.3728 Jinke Rd,
Pudong New District
Shanghai, 201203- China
Tel.: +400 821 8006
Email: info_china@vwr.com

Czech Republic

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Strábrná Skalce
Tel.: +420 321 570 321
Email: info.cz@vwr.com

Denmark

VWR International A/S
Tobaksvejen 21
2860 Søborg
Tel.: +45 43 86 87 88
Email: info.dk@vwr.com

Finland

VWR International Oy
Valimotie 17-19
00380 Helsinki
Tel.: +358 9 80 45 51
Email: info.fi@vwr.com

France

VWR International S.A.S.
Immeuble Estréo
1-3 rue d'Aurion
93114 Rosny-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00
(international)
Email: info.fr@vwr.com
* 0,18 € TTC/ min + prix appel

Germany

VWR International GmbH
Hilpertstrasse 20a
D - 64295 Darmstadt
Freecall: 0800 702 00 07
Tel.: +49 (0) 6151 3972 0 (international)
Email: info.de@vwr.com

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 (52) 521-130
Email: info.hu@vwr.com

India

VWR Lab Products Private Limited
No.139. BDA Industrial Suburb,
6th Main, Tumkur Road,
Peenya Post,
Bangalore, India – 560058
Tel.: +91-80-28078400
Email: vwr_india@vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR
International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 01 88 22 222
Email sales.ie@vwr.com

Italy

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02-3320311/02-487791
Email: info.it@vwr.com

Korea

VWR International ~
17, Daehak 4-ro, Yeongtong-gu,
Suwon-si, Gyeonggi-do
Tel.: +82 31 645 7256
saleskorea@avantorsciences.com

Middle East & Africa

VWR International FZ-LLC
DSP Laboratory Complex
125, Floor 01
Dubai, United Arab Emirates
Tel.: +971 4 5573271
info.mea@vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 020 4808 400
Email: info.nl@vwr.com

Norway

VWR International AS
Brynsalleen 4
0667 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 058 32 38 210
Email: info.pl@vwr.com

Portugal

VWR International - Material de
Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, n° 6
2610-088 Amadora
Tel.: +351 21 3600 770
Email: info.pt@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Email: sales.sg@vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
Email: info.es@vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 08 621 34 00
Email: info.se@vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 044 745 13 13
Email: info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
Email: uksales@vwr.com