IKA® RV 10 control

Betriebsanleitung	DE	5
Operating instructions	EN	38
Mode d'emploi	FR	71
Instrucciones de manejo	ES	104
Инструкция по эксплуатации	RU	137

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Contents

	Page
EC-Declaration of conformity	3
Explication of warning symbols	38
Warranty	38
Safety instructions	39
Correct use	40
Unpacking	41
Useful information	41
Vacuum control: two-position and speed-control	42
Setting up	43
Drive RV 10 control	
Heating bath	
Glassware	
Hose system	
Interfaces and outputs	51
Commissioning	53
Maintenance and cleaning	65
Accessories	65
Error codes	68
Technical data	69
Solvent table (excerpt)	70

Explication of warning symbols

General hazard

Danger of being burnt!

Warranty

In accordance with **IKA**[®] warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Safety instructions

For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations. **When working under a vacuum in particular!**
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - splashing liquids,
 - body parts, hair, clothing and jewellery getting caught,
 - injury as a result of glass breakage.
- **Caution!** Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
 - stress due to incorrect assembly,
 - external mechanical hazards,
 - local temperature peaks.
- Ensure that the stand does not start to move due to vibrations respectively unbalance.

- Beware of hazards due to:
 - flammable materials,
 - combustible media with a low boiling temperature.

- **Caution!** Only process and heat up media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set. The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used.
- Do **not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Tasks with the device must only be performed when operation is monitored.
- Operation with excess pressure is not permitted (for cooling water pressure see "Technical Data")
- Do not cover the ventilation slots of the device in order to ensure adequate cooling of the drive.
- There may be electrostatic discharges between the medium and the drive which could pose a direct danger.
- The appliance is not suitable for manual operation.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Refer to the operating instructions for the HB 10 heating bath.
- Refer to the operating instructions for the accessories, e.g. vacuum pump.
- Position the positive pressure outlet of the vacuum pump under a fume hood.
- Use the HB 10.2 cover or the HB 10.1 spray guard.

- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Adapt the quantity and the type of distilland to the size of the distillation equipment. The cooler must work properly. Monitor the coolant flow rate at the cooler outlet.
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at cooler) in order to prevent a pressure build-up.
- Please note that dangerous concentrations of gases, vapours or particulate matter can escape through the outlet at the cooler. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.
- Evacuated glass vessels must not be heated on one side; the evaporating flask must rotate during the heating phase.
- The glassware is designed for operation under a vacuum of up to 1 mbar. The equipment must be evacuated prior to heating (see chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, uncondensed vapours must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief
- **Caution!** Avoid peroxide formation. Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing! Keep liquids that tend to form organic peroxides away from light, in particular from UV rays and check them prior to distillation and exhaust for the presence of peroxides. Any existing peroxides must be eliminated. Many organic compounds are prone to the formation of peroxides e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.
- **Danger of burning!** The heating bath, tempering medium, evaporator piston and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.
- **ATTENTION!** Avoid delayed boiling! Never heat the evaporating flask in the heating bath without switching on the rotary drive! Sudden foaming or exhaust gases indicate that flask content is beginning to decompose. Switch off heating immediately. Use the lifting mechanism to lift the evaporator piston out of the heating bath. Evacuate the danger zone and warn those in the surrounding area!
- **Caution!** Never operate the device when the evaporator piston is rotating and the lift is raised. Start at low speed and lower the evaporating flask into the heating bath, increase to desired speed. Otherwise hot tempering medium may be sprayed out!
- Set the speed of the drive so no tempering medium is sprayed out as a result of the evaporator piston rotating in the heating bath. If necessary reduce the speed.
- Do not touch rotating parts during operation.
- Imbalance may result in uncontrolled resonance behaviour of the device or assembly. Glass apparatus may be damaged or destroyed. In the event of unbalance or unusual noises, switch off the appliance immediately or reduce the speed.
- A vacuum may form inside the glassware in case of power outage. The device will vent automatically when it is switched on again.
- The appliance does not start up again automatically following a cut in the power supply.
- The device is only disconnected from the power supply network if the device power switch is off or the plug is pulled out.
- The socket for the mains cord must be easily accessible.

- Automatic operation: Always check evaporation processes before starting automatic operation. Automatic operation must not be used for unfamiliar evaporation processes.

Safety lift

When the device is switched off or the power supply disconnected, the internal safety lift removes the evaporating flask from the heating bath.

When the safety lift is powered off, it is designed for a maximum total weight (glassware and solvent) of 3.1 kg.

Sample calculation of total weight with vertical glassware and a 1 litre flask:

Condenser + receiving flask + evaporating flask + fittings =
1200 g + 400 g + 280 g + 100 g = 1980 g

Maximum solvent load = 3100 g – 1980 g = 1320 g

Due to the design, the safety lift cannot be guaranteed for higher loads!

When using other types of condensers such as dry ice or intensive condensers as well as when using return distillation distributors with slip-on condensers, it may be necessary to reduce the load by the added weight of the glass apparatus.

Thus, prior to distillation, check whether the lift goes up without power when laden with the glassware and distillation material.

The safety lift must be checked daily prior to operation. Use the motor to move the lift to the lowest position and press the "Power" key on the

For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- The appliance may only be opened by experts.

The following automatic modes are available in addition to manual mode:

- 100% distillation,
- Volume-/quantity-dependent distillation.

front plate or the main switch on the back right side of the device using the maximum total weight of 3.1 kg.

➔ The evaporating flask is lifted out of the heating bath.

When using for the first time or after a long period of inactivity, proceed as follows: use the motor to move the lift several times to the top and bottom positions by pressing the lift function keys prior to interrupting the power supply.

Note: Switching off and disconnecting the power supply shows a delay compared to switching off using the "Power" switch on the front panel.

If the safety lift is not working, please contact the **IKA®** Service department.

For the evaporating equipment (evaporating flask plus contents), the maximum permissible weight is 3.0 kg! With larger loads there is the risk of broken glass on the steam pipe!

Ensure that the safety lift has been powered down.

When working with large loads, always use low speeds. Unbalanced loads can result in a broken steam pipe!

- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids etc.
- Protect the appliance and accessories from bumps and impacts.

Correct use

• Use

Together with the accessories recommended by **IKA®**, the device is suitable for:

- quick and gentle distillation of liquids
- evaporation of solutions and suspensions
- crystallisation, synthesis or cleaning of fine chemicals
- drying of powder and granulate material
- recycling of solvents

Mode of operation: Tabletop device

• Range of use

- | | |
|----------------|----------------|
| - Laboratories | - Schools |
| - Pharmacies | - Universities |

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications or if the appliance or the printed circuit board are modified by third parties.

Unpacking

• Unpacking

- Please unpack the device carefully
- In the case of any damage a fact report must be sent immediately (post, rail or forwarder)

• Delivery scope

see table

	Drive RV 10 control	Heating bath HB 10	Vertical glassware RV 10.1	Vertical glassware RV 10.10 coated	Woulff bottle	Bracket	Cooler locking device, compl.	Vacuum hose (2 x 0.55 m)	Water discharge hose (1 x 1 m)	Operating instructions	Ring spanner	Handhold	RV 10.4002 vacuum valve for laboratory vacuum system	N 920 speed-controlled vacuum pump
RV 10 control V	x	x	x		x	x	x	x	x	x	x	x	x	
RV 10 control V-C	x	x		x	x	x	x	x	x	x	x	x	x	
RV 10 control FLEX	x	x			x	x	x	x	x	x	x	x	x	
RV 10 control V auto	x	x	x		x	x	x	x	x	x	x	x		x
RV 10 control V-C auto	x	x		x	x	x	x	x	x	x	x	x		x
RV 10 control FLEX auto	x	x			x	x	x	x	x	x	x	x		x

Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.

The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C). Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation cooler should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a Woulff bottle and/or a vacuum separator. Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems.

Speed, temperature, piston size and system pressure all affect the evaporator air capacity.

The optimum capacity of the flow-through condenser is approx. 60%.

This corresponds to condensation on approx. 2/3 of the cooling coil. With larger capacities there is the risk that the uncondensed solvent vapour will be extracted.

In order to prevent this, the condenser loading safety monitoring system can be activated. Before activating the safety monitoring, first carry out several distillations, noting the value given under menu "Settings", menu option "Drive", item "Maximum heat rejection rate". Now calculate the heat rejection rate using the following formula:

Heat rejection rate threshold = maximum heat rejection rate * 1.3

Enter this value under menu "Settings", menu option "Distillation", item "Heat rejection rate threshold". The device will now stop the measurement and issue an error message if the threshold for the condenser is exceeded. If the configuration of the condenser equipment or the distillation process is modified, then it may be necessary to re-evaluate the maximum condenser loading.

The factory setting for this value is 900 W.

The device is equipped with a piston-operated safety mechanism. The glass apparatus may contain a vacuum following interruption to the power supply; vent the system before switching back on. If the power cuts out, the evaporator piston is automatically lifted out of the heating bath by an integrated gas spring.

CAUTION! The safety lift must be checked daily prior to operation. For more information see the section entitled "Safety Instructions - Safety Lift"!

In addition to offering a full range of manual and semi-automatic evaporation operating modes, the RV 10 control rotary evaporator can also be used for fully automatic and quantity-controlled evaporation processes. Standard fittings for fully-automatic operation include a vacuum controller, a cooling water differential temperature measurement sensor and a cooling water flow rate meter. The vacuum can be controlled in both two-position and speed-controlled pumping operation with the help of the integrated vacuum controller. The device is designed for operation with a cooling water supply system (e.g. laboratory thermostat), but can also be run off a water supply line. Please refer to the Technical Data for information on cooling water pressure, temperature constancy, flow rate and local regulations.

We recommend using the RV 10.5001 water regulator valve when operating directly off a water supply line.

This valve allows the cooling water flow rate to be adjusted and automatically cuts off the cooling water supply on completion of the distillation process.

The RV 10.5002 water filter can be used to prevent particles of dirt passing from the water supply line into the valve mechanism.

If pressure reduction is required, we recommend fitting the RV 10.5003 pressure regulator valve in the pipe immediately downstream of the pick-off point.

Automatic volume-dependent distillation: The apparatus must be brought up to operating temperature. This is achieved by conducting a test distillation.

Vacuum control: two-position and speed control

The RV 10 control rotary evaporator allows you to set the desired vacuum using two-position and speed control with the help of the integrated vacuum controller and pressure sensor.

Two-position control

A vacuum is made in the glassware with the help of a vacuum pump. The vacuum pump operates at constant speed which generally cannot be adjusted. Once the target value setting is reached, the RV 10.4002 vacuum valve interrupts the suction line (included in scope of supply). If additional switching on and off of the vacuum pump is desired, please use the RV 10.4003 pump control (accessory).

The device will fall slightly short of the target value setting due to the time difference in the millisecond range with the vacuum value, the target value comparison and the switching of the vacuum valve. Natural leakage from the apparatus now causes another increase in the pressure value of the system, which in turn is detected by the integrated pressure sensor.


If the target value is exceeded, the vacuum valve is reopened and the vacuum pump in operation lowers the pressure. The pressure curve shown on the display fluctuates around the set target value.

The pressure difference between activating/deactivating the valve can be set using the vacuum hysteresis value (see Vacuum Settings). The value is set to 10 mbar as standard.

The quality of the control depends on the suction power of the pump (rotation speed), the tightness of the system and signal delay time of the electronics used.

Automatic boiling point recognition is not possible when using two-position control.

Schematic view of a two-position vacuum control

Speed-vacuum control

Imprecision of the two-position control is avoided by using the speed-vacuum control.

However a vacuum pump with a controlled speed is needed for speed-vacuum control; the RV10.4002 vacuum valve is not required and must not be connected!

To do so, connect the vacuum pump (e.g. KNF N920 speed controlled pump) to the port on the rear of the rotary evaporator.

With this kind of control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value.

Once the target value is reached, the pump only operates according to the leakage rate.


It is possible to achieve quieter operation and more exact vacuum control.

The speed-vacuum control on the RV 10 is controlled automatically as soon as an appropriate vacuum pump is attached.

Automatic boiling point recognition is possible with this type of control, i.e. the system reaches and keeps the solvent at its boiling point in automatic mode.

As it is not necessary to know the boiling point of the solvent when in automatic boiling point recognition mode (unlike in the case of volume-dependent distillation), these two operating modes are mutually exclusive.

Schematic view of speed-vacuum control

Setting up

Drive RV 10 control

Caution! Loosen transportation lock (Fig. 4a)


Fig. 4a

- ☞ Hold the lift with your hand and remove the thumb screw on the back of the appliance (A).
- ☞ Once the transportation lock has been removed, the lift moves slowly to its upper end position. The distance is approx. 140 mm.
- ☞ Connect the device to the power supply (B) using the supplied power cord (B).

Adjustable base (Fig. 4b)


Fig. 4b

Mounting the Woulff bottle (Fig. 5)

- ☞ Fit the hose connecting piece (C) on the left side of the lift.
- ☞ Mount the holder (D) onto the hose connecting piece (C).


Fig. 5

- ☞ Insert the bottle and attach the supplied hose connectors to the bottle (Fig. 6).

Fig. 6

- ☞ Remove the clamping device for the angle setting of the rotation drive on the right side of the lift by rotating the knurled screw counter clockwise (by lightly pressing and turning at the same time, the knurled screw comes out farther).
- ☞ Set the drive at an angle of approx. 30° (Fig. 7).


Fig. 7

- ☞ Then secure the rotation drive from being accidentally turned by tightening the knurled screw in a clockwise direction.

Heating bath HB 10

Caution! Refer to the chapter entitled "Commissioning" in the heating bath instruction manual!

- ☞ Place the heating bath on the stand of the rotation drive and push it into the left position (Fig. 8).


Fig. 8

Note: Data is exchanged between the drive unit and the heating bath by means of an infrared link (F). Please note that reliable communication is only guaranteed when the infrared beam has clear line-of-sight to the detector.

Glassware

- ☞ Open the locking device of the steam pipe by turning it 60° counter clockwise (D), (fig. 9).


Fig. 9

- ☞ Feed the steam pipe in until it stops.
- ☞ Then lock the locking device by turning it clockwise by 60° (Fig. 10).

- ☞ Check the axial locking device on the steam pipe.

Push-off mechanism

The plastic screw nut (H) helps loosen tight-fitting piston ground-in connections (Fig. 10).

- ☞ Hold the locking device (G) and loosen the plastic screw nut (H).
- ☞ To do so, hold the tight-fitting evaporator piston and turn the plastic screw nut (H) until the evaporator piston neck.

Note: Prior to commissioning, hand-tighten the plastic screw nut (H) left-aligned. This will hold the steam pipe locking device (G) firmly in place.

Fitting the condenser seal

- ☞ Place the RV 10.8001 condenser seal in the condenser receptacle and fit the glassware to the device according to the assembly instructions (Fig. 11 a,b, c and d).

Fig. 11a

Fig. 11b

First use - Fitting the seal RV 10.8001

- ☞ Insert vapour tube (1).
- ☞ Insert seal RV 10.8001 (2).
- ☞ Slide the union nut (3b) over the flange on the condenser (3a).
- ☞ Also slide the annular spring (3c) over the flange on the condenser (3a).
- ☞ Position the condenser (3a) on top of the seal (2).
- ☞ Tighten the cap nut hand-screwed (3b).

Note: Follow the mounting instructions for the glassware


Fig. 11c

Assembling the vertical glassware cooler locking device

- ☞ Assemble the condenser locking device according to the diagram (fig. 12).
- ☞ Mount the plate (I) using the thumb screw (J).
- ☞ Put the support rod (L) on the plate (I) and attach it with the screw nut (K).
- ☞ Attach the rubber protector (M).
- ☞ Fasten the Velcro® band (N) to the support rod (L).
- ☞ Secure the vertical glassware with the Velcro® band (N).



Fig. 12

Glassware

Mounting the glassware

Item	Designation	Quantity	
		diagonal glassware	vertical Glassware
1	Receiving flask	1	1
2	Clamp NS 29 (stainless steel)	1	1
3	Steam pipe	1	1
4	Stopcock	1	1
5	Pipe tube	1	1
6	Condenser	1 Diagonal condenser	1 Vertical condenser
7	Connection	1 Introduction sleeve	1 Vacuum connection
8	Clamp NS 29 (plastic)	1	1
9	Seal RV 10.8001	1	1
10	Evaporator piston 1000 ml	1	1
11	Ball joint clamp RV 05.10	1	1
12	Screw joint cap	4	4
13	Hose connection	4	4
15	Locking nut	1	1
16	Spring ring	1	1

Removing the condenser

- ☞ Use the ring spanner provided to loosen union nuts that are tightly fitted.
- ☞ Loosen the union nut by turning anticlockwise.
- ☞ Remove the Velcro®.

Fig. 13

Description of special condensers

- **RV 10.3 Vertical-intensive condenser with manifold**
This vertical-intensive condenser features a double jacket design for particularly efficient condensation.
Also available with coating (RV 10.30)
- **RV 10.4 Dry ice condenser**
Dry ice condenser for distilling low-boiling solvents.
Cooling by dry ice, no cooling water required. Maximum condensation thanks to low temperatures.
Also available with coating (RV 10.40)
(Cannot be used in automatic mode for the RV10 control)
- **RV 10.5 Vertical-condenser with manifold and cut-off valve for reflux distillation**
Also available with coating (RV 10.50)
- **RV 10.6 Vertical-intensive condenser with manifold and cut-off valve for reflux distillation**
This vertical-intensive action condenser features a double jacket design for particularly efficient reflux distillation.
Also available with coating (RV 10.60)

Hose system

Schematic view of hose system (rear view)


Fig. 14

Vacuum - two-position control

☞ Fit one of the optional RV 10.4001, RV 10.4002 or RV1 0.4003 vacuum valves, as appropriate for the type of vacuum system used, in the bracket provided for this purpose and connect the vacuum hose to the valve (Figure 15a, b, c).

Connection RV 10.4001



Fig. 15a

Connection RV 10.4002



Fig. 15b

Connection RV 10.4003



Fig. 15c

Vacuum - speed-controlled operation

☞ No additional valve is required for vacuum normal operation with a speed-controlled pump. Connect the pump's vacuum hose directly to the pressure sensor hose connection. Connect the electrical RV10 control cable to the controller (Fig. 15e).

Fig. 15d

Fig. 15e

Connection RV 10 temperature sensor (dT)

Connect the temperature sensor to the lower socket (dT). Always check for correct direction of the arrow.


Fig. 15f

dT

Water

☞ Connect the water inlet to the water supply (Figure 16a; Figure 16b with optional RV 10.5001 water regulator valve). Please refer to the technical specifications relating to the water supply.

The RV 10.5001 water regulator valve is not suitable for use on the condenser unit because it would reduce the flow rate too greatly


Fig. 16a


Fig. 16b

☞ Connect the water hoses to the glass condenser (short hose (O) = lower outlet, long hose (P) = upper inlet) and fix hoses in place (fig. 17).

Note: For special condensers see "Description of special condensers" section.


Fig. 17

☞ Connect the water outlet hose (supplied) by pushing it onto the nipple until the stop inside the connection socket is reached (fig. 18).

Caution! Ensure that the condenser inlet and condenser outlet are connected correctly.

Automatic distillation cannot be carried out if the water inlet and outlet hoses are not fitted according to the instructions and illustrations acc. fig. 14, since this would cause the temperature readings for the inlet and outlet to be measured incorrectly.


Fig. 18

☞ Release the hose connection using the supplied tool.

Fig. 19

Schematic view of connections (rear view)

☞ Insert the valve connector cable (RV 10.5001 or RV 10.4001/2 resp. RV 10.4003 Valve) resp. the network cable (RV 10.4003 pump control) into the appropriate socket (fig. 20) or connect the speed-regulated vacuum pump. As a result, the RV10 control automatically switches to normal speed-vacuum control operation.


Fig. 20

☞ Connect the pressure sensor to the Woulff bottle and the Woulff bottle to the vacuum connector on the condenser using the vacuum hoses supplied. Please note that the vacuum must always be connected to the highest connection point on the condenser (fig. 21).


Fig. 21

HB 10

Data is exchanged between the heating bath and the drive unit by means of an infrared link. These are located on the left display side of the heating bath or on the right side of the drive unit. Do not place any objects between the two operating units as otherwise the data transfer may be interrupted (fig. 22)!

Fig. 22

RV10 control

The device can be operated in "Remote" mode via an RS 232 or USB interface using labworldsoft® laboratory software.

The RS 232 interface at the back of the device, see fig. 20, is fitted with a 9-pole SUB-D jack and can be connected to a PC. The pins have serial signals. The USB interface is located on the left-hand side of the display on the drive unit and can be connected to a PC using the USB cable supplied.

Note: Please note the system requirements as well as the operating instructions and help section included with the software.

USB Interface

The Universal Serial Bus (USB) is a serial bus system used to connect the RV 10 control with the PC (Fig. 23). Devices that support USB can be connected to each other whilst they are running (hot pluggable) and provide automatic recognition of the connected devices and their properties.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode to update the solvent library. Navigate to <http://www.ika.net/ika/lws/download/RV10.cfg> to update.

Fig. 23

Installation

Connect the RV 10 control to the PC using the USB data cable. The RV 10 control will then transmit information to the Windows operating system to tell it which device drivers are required. Windows will then either:

- Load the driver,
- Install the driver automatically if not already installed,
- Prompt the user to perform a manual installation.

Navigate to <http://www.ika.net/ika/lws/download/stmcddc.inf>.

Serial interface RS 232 (V24)

Configuration

- The functions of the interface lines between the device and the automation system are a selection from the signals specified in the EIA standard RS232 C, as per DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 C applies in accordance with DIN 66 259 part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s.
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request.

Command syntax and format

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the device (Slave).
- The Rotary evaporator sends only at the computer's request. Even fault indications cannot be sent spontaneously from the Rotary evaporator to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a number is a dot (Code: hex 0X20E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific IKA® commands serve only as low level commands for communication between the rotary evaporator and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the rotary evaporator. The IKA® software package, labworldsoft, provides a convenient tool for controlling rotary evaporator and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarises the (NAMUR) commands understood by the IKA® control equipment.

Abbreviations used:

m =	Numbering parameter (integer)
X = 2	Temperature heating bath
X = 3	Safety temperature heating bath
X = 4	Speed
X = 60	Interval time (1-99 seconds, 1 <= m >=99)
X = 61	Timer (1-199 minutes, 1 <= m >=199)
X = 62	Upper lift position (OUT_SP_62 1-> drive lift up)
X = 63	Lower lift position (OUT_SP_62 1-> drive lift down)
X = 66	Value vacuum controller
X = 70	Vacuum controller hysteresis
X = 74	Tempering medium (OUT_SP_74 0=Öl, OUT_SP_74 1=water)

NAMUR Commands	Function
IN_NAME	Request designation
IN_PV_X X = 4	Read actual value
IN_SOFTWARE	Request software Id-number, date and version
IN_SP_X X = 4	Read target value input
OUT_SP_X m X = 1,60,61,62	Set target value to m
RESET	Switch to normal operation
START_X X = 4,60,61,62	Switch on appliance (remote) function
STATUS	Status output 0:Manual operation without interruption 1:Automatic operation Start (without interruption) ERROR z (z error number see table)
STOP_X X = 1,60, 61, 62	Switch off appliance function. The variable set with OUT_SP_X remains

PC 1.1 Cable (Fig. 24)

This cable is required to connect the 9-pin connector to a PC. (fig. 24).

Fig. 24

Commissioning

Basic guide to using the menu system

Selecting a menu

- ☞ Turn the Rotating/ Pressing knob to the right/left to select the desired menu option.

Note: The active menu item is shown with a black background on the display.

Editing values

- ☞ Press the Rotating/ Pressing knob
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to complete the process.

Exiting a menu

- ☞ Press the SET button.
- ☞ Current values are stored.
- Next-highest menu level is displayed.

- ☞ Press the ESC button.
- ☞ Current values are not stored.
- Next-highest menu level is displayed.

Note: Some menu items do not contain values that can be stored. Select SET or ESC to exit the menu.

Note: Certain menu items have functionality that varies from or adds to that described above.

Display screen during the distillation process


Each operating mode has its own dedicated display screen. The following display characteristics are common:

- When a distillation is in progress it is not possible to leave the active display indicated.
- Press SET or ESC to exit the active display mode when no distillation is in progress. This will return the display to the main menu.

Error messages

- ☞ Press ESC to acknowledge the fault if an error message appears on the display.
- ☞ The error message will disappear.

Note: In the case of serious faults, an additional display message will appear indicating that the device must be switched off and only switched back on again when the fault has been rectified.

The unit is ready for service when the mains plug has been plugged in.

Switching on the device

- ☞ Switch on the device on the right side of the device by the mains plug (fig.25).
- ☞ Device functions activated.

Fig. 25

Note: The heating bath must also be switched on.

Functional description (in as-delivered state)


Factory settings

The values specified below correspond to the as-delivered state (display output/factory language setting is English)


- ☞ The system check is performed while the start screen is displayed; this takes no more than 30 seconds.

- ☞ The main menu will appear after a few seconds.

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Confirm the selection by pressing the Rotating/ Pressing knob.

"Settings" menu

The following options are available from the "Settings" menu.

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Confirm the selection by pressing the Rotating/ Pressing knob.

Distillation settings

Settings Distillation	
Type of distillation	100%
Drying	Off
Limit of cooler power	900 W
Maximum cooler power	0 W
Flow Check	On
Unit Amount	mL
Boiling Detection	No
SET = Save ESC = Cancel	

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to complete the process.

SET key: Exit menu and save changes.

ESC key: Exit menu without saving changes.

Distillation method

Select between the automatic distillation methods "Volume" and "100%". This selection only has an effect during auto-distillation.

Volume

The measured values of the cooling water flow rate and the cooling water temperature difference form the basis for a heat balance calculation used to determine the quantity of distillate at every stage of the distillation. The distillation process is halted when the specified quantity of distillate has been reached.

100%

The distillation process is stopped when the measured cooling water temperature difference falls below a threshold value, i.e. as soon as a solvent has been completely distilled off.

Drying

The cooling water temperature difference is not monitored if this option is activated, e.g. use in processes for drying powdered media.

Heat rejection rate threshold

A limit value can be set for the heat rejection rate that is appropriate for the condenser used; see "useful information" section. The actual heat rejection rate is calculated during every distillation process. If the threshold is exceeded, the distillation will be stopped and an error message displayed.

Maximum heat rejection rate

The maximum heat rejection rate achieved during the last distillation is displayed. This value is for information purposes only.

Quantity meas. unit

Select between millilitre and gram as the unit of quantity.

Boiling point recognition

Boiling point recognition can be either activated or deactivated. Boiling point recognition is only possible in "100%" mode for the distillation method.

The boiling point recognition function automatically detects the boiling point of a solvent.

Vacuum settings

Settings Vacuum	
Hysteresis	10 mbar
Scaling time	0 min
Scaling vacuum	0 mbar
Unit Vacuum	mbar
SET = Save ESC = Cancel	

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to complete the process.

SET key: Exit menu and save changes.

ESC key: Exit menu without saving changes.

Hysteresis

The default settings for the hysteresis (as shown) are suitable for use in most situations. The (vacuum) hysteresis value describes the difference between the closing and opening pressures for the vacuum valve.

Time scaling

The time axis of the vacuum-rotation-time graph is scaled using the value specified here. If the value = 0, then the time axis is scaled automatically.

Vacuum scaling

The vacuum axis of the vacuum-rotation-time graph is scaled using the value specified here. If the value = 0, then the vacuum axis is scaled automatically.

Vacuum meas. unit

Select between the units of measurement for a vacuum: mBar, Torr, and hPascal.

Drive settings

Settings Drive	
Right/Left interval	0 sec
Timer	0 min
Scaling time	0 min
Scaling rotation	0 rpm
SET = Save	ESC = Cancel

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to complete the process.

SET key: Exit menu and save changes.

ESC key: Exit menu without saving changes.

Right/left interval

Setting a particular value X for the right/left interval will cause the drive to change the direction of rotation every X seconds.

Note: In interval mode the maximum speed is limited to 200 rpm.

Timer

The timer value determines the time after which a manual distillation will be stopped.

Time scaling

The time axis of the vacuum-rotation-time graph is scaled using the value specified here. If the value = 0, then the time axis is scaled automatically.

Rotation scaling

The rotation axis of the vacuum-rotation-time graph is scaled using the value specified here. If the value = 0, then the rotation axis is scaled automatically.

Bath settings

Settings Heating Bath	
Bath medium	Water
SET = Save	ESC = Cancel

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to complete the process.

SET key: Exit menu and save changes.

ESC key: Exit menu without saving changes.

Medium bath

Select either water or oil as the medium for the heating bath.

Language

Language	
English	
Deutsch	
Français	
Español	
中文	
日本語	
Polski	
SET = Save	ESC = Cancel

- ☞ Turn the Rotating/ Pressing knob to the left/right to select the language.
- ☞ Confirm the selection by pressing SET.

Service

Service			
RV10.5001	Off	38.0L/h	
RV10.4001/2	Off		
RV10.4003 Valve	Off		
V4 Vacuum venting	Off		
Pump	Off		
Adjustment	3.2 K		
Operation hours	14 h		
ESC = Back			

- ☞ Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to complete the process.

ESC key: Exit menu, all changes to switch settings are reset, the calibration remains unaffected.

Switching valves and pump

To check functionality for servicing purposes, use the Service menu to activate/deactivate directly the valves included in delivery or optional and the pump, e.g. RV 10.5001.

The RV 10.5001 valve switches the water circuit on and off. This row also shows the current cooling water flow rate.

Calibrate

The "Calibrate" option is used for calibration of the temperature sensor.

The temperature difference displayed will be calibrated to zero. This action cannot be reversed.

The temperature sensor must be calibrated whenever it is replaced or changed, or if a new temperature sensor is installed.

Please contact the Service department.

The temperature sensor supplied with the system is calibrated ex-works.

Performing calibration

- ☞ Activate the "Calibrate" menu option by turning the Rotating/ Pressing knob.
- ☞ Remove both sensors. To remove a sensor, push back the outer ring to release the locking device on the push-in connector, while at the same time pulling the temperature sensor out of the connector, see fig. 26.

Fig. 26

- ☞ To reconnect a temperature sensor, push it into the push-in connector until the limit stop is reached; an initial resistance will be felt.

- ☞ Fill water at room ambient temperature into a glass beaker (approx. 500 ml). Fully submerge both sensors into the water, see fig. 27.

Fig. 27

- ☞ Wait until the temperature display in the "Service" menu, menu option "Calibrate" has stabilised, e.g. 0.2K

- ☞ Press the Rotating/ Pressing knob.
⇒ The temperature difference displayed will be calibrated to zero.


Calibration is only possible within the temperature range +0.5 K to -0.5 K. If the temperature is outside this range, then either a temperature sensor is faulty or connected incorrectly. Please contact the Service department.

If the characteristics of the cooling water supply deviate significantly from the specifications given in the "Technical specifications", with the result that the initial conditions for automatic distillation are not achieved, then it may be necessary to recalibrate the temperature sensors. This may, however, lead to a reduction in the accuracy with which the distillate quantity is measured in "Auto-distillation" mode.

Operation hours

This value serves as information and cannot be changed.

"Manual mode without boiling point recognition" menu

"Manual mode without boiling point recognition" display screen

1. Operating mode display
2. Vacuum display (actual value)
3. Vacuum display (target value)
4. Control field for „Start/ Stop“ distillation
5. Timer display
6. Control field for storing the distillation sequence
7. Vacuum-rotation-time graph
8. Display of flow rate, differential temperature, and condenser heat rejection rate
9. Control field for rotary speed display (target value) and speed measurement unit
10. Error message in case of fault, otherwise status indicator

Setting the rotary speed

- ☞ Press the Rotating/ Pressing knob to change the target value.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to save and complete the process.

Note: When you select the speed > 100 rpm, smooth start is automatically activated.

Setting the vacuum

- ☞ Press the Rotating/ Pressing knob to change the target value.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to save and complete the process.

Start

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Start".
- ☞ Press the Rotating/ Pressing knob to start the distillation.
- ☞ The display field will change to "Stop".


Note: If the timer has been configured, then the time will start to count backwards. If a timer has not been configured then the time elapsed since the start is shown.

Stop

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Stop".
- ☞ Press the Rotating/ Pressing knob to stop the distillation.
- ☞ The display field will change to "Start".

Saving the distillation sequence as a procedure on completion of distillation


- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Save".
- ☞ Press the Rotating/ Pressing knob to display the "Procedures" screen.

- ☞ Turn the Rotating/ Pressing knob to the left/right to select one of the ten procedures. The sequence will be shown graphically.
- ☞ Press SET, to overwrite the selected procedure with the current procedure (the procedure that has just been completed) and exit the menu. Pressing ESC will prevent the selected procedure from being overwritten.

Note: Distillation sequences saved in this manner can be repeated at a later date in the user-defined mode by selecting the corresponding procedure.

"Manual mode with boiling point recognition" menu

"Manual mode with boiling point recognition" display screen

1. Operating mode display
2. Vacuum display (actual value)
3. Vacuum display (target value)
4. Control field for „Start/ Stop“ distillation
5. Timer display
6. Control field for storing the distillation sequence
7. Vacuum-rotation-time graph
8. Display of flow rate, differential temperature, and condenser heat rejection rate
9. Control field for rotary speed display (target value) and speed measurement unit
10. Error message in case of fault, otherwise status indicator

Setting the rotary speed

- ☞ Press the Rotating/ Pressing knob to change the target value.
- ☞ Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned).
- ☞ Press the Rotating/ Pressing knob again to save and complete the process.

Note: When you select the speed > 100 rpm, smooth start is automatically activated.

Start

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Start".
- ☞ Press the Rotating/ Pressing knob to start the distillation.
- ☞ The display field will change to "Stop".


Note: If the timer has been configured, then the time will start to count backwards. If a timer has not been configured then the time elapsed since the start is shown.

Stop

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Stop".
- ☞ Press the Rotating/ Pressing knob to stop the distillation.
- ☞ The display field will change to "Start".

Saving the distillation sequence as a procedure on completion of distillation


- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Save".
- ☞ Press the Rotating/ Pressing knob to display the "Procedures" screen.

- ☞ Turn the Rotating/ Pressing knob to the left/right to select one of the ten procedures. The sequence will be shown graphically.
- ☞ Press SET, to overwrite the selected procedure with the current procedure (the procedure that has just been completed) and exit the menu.
Pressing ESC will prevent the selected procedure from being overwritten.

Note: Distillation sequences saved in this manner can be repeated at a later date in the user-defined mode by selecting the corresponding procedure.

"Auto-distillation mode without boiling point recognition" menu


"Auto-distillation mode without boiling point recognition" display screen

1. Operating mode display
2. Control field for solvent
3. Vacuum display (actual value)
4. Vacuum display (target value)
5. Control field for „Start/ Stop“ distillation
6. Timer display
7. Distillation progress indicator (% of required quantity of distillate)
8. Vacuum-rotation-time graph
9. Display of flow rate, differential temperature, and condenser heat rejection rate.
10. Rotary speed display (target value) and speed measurement unit
11. Error message in case of fault, otherwise status indicator

Changing the selected solvent


- ☞ Turn the Rotating/ Pressing knob to the right/left to select the control field "Solvent".
- ☞ Press the Rotating/ Pressing knob to display a list of the various solvents available.


- ☞ Turn the Rotating/ Pressing knob to the right/left to select a "Solvent".
- ☞ Press the SET key.
- ☞ The solvent selected will be used for the subsequent distillation.

Note: The display shows the parameters for each of the solvents. The parameters can only be changed for the user-definable solvents UserSolvent1 ... UserSolvent5. Please refer to the next section for instructions on how to do this.

User defined solvents

- ☞ Select one of the options NewSolvent1 ... NewSolvent5 from the "Solvents" screen.
- ☞ Select one of these solvents and press the Rotating/ Pressing knob.


All the parameters displayed for the solvent can now be modified.

- ☞ Turn the Rotating/ Pressing knob to the right/left to select the parameter to be changed.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to the right/left to change the selected parameter.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Press SET to save all values.
- ☞ This action selects NewSolvent1 as the current solvent.

Note: It is not possible to return to the selection list.

"Auto-distillation mode with boiling point recognition" menu



"Auto-distillation mode with boiling point recognition" display screen

1. Operating mode display
2. Control field for solvent
3. Vacuum display (actual value)
4. Vacuum display (target value)
5. Control field for „Start/ Stop“ distillation
6. Timer display
7. Distillation progress indicator (% of required quantity of distillate)
8. Vacuum-rotation-time graph
9. Display of flow rate, differential temperature, and condenser heat rejection rate
10. Rotary speed display (target value) and speed measurement unit
11. Error message in case of fault, otherwise status indicator

Note: volume distillation is not possible as the solvent must be known in this mode.


This mode allows selection of the bath temperature and speed of rotation via the "Boiling point recognition" control field instead of selection of the solvent.

Starting automatic distillation

Press the Rotating/ Pressing knob with the "Start" or "Continue" field selected. If the "100% distillation" or "Drying" options are active, then the distillation will start immediately, provided that the necessary conditions are fulfilled.

When the "Volume" distillation method is used, a further screen with additional distillation parameters will be displayed.

- Turn the Rotating/ Pressing knob to the right/left to select a menu option.
- Press the Rotating/ Pressing knob.
- Turn the Rotating/ Pressing knob to change the value (the magnitude of the change is dependent on the speed at which the knob is turned)
- Press the Rotating/ Pressing knob again to complete the process.

Distillation

This field indicates the distillation method selected in Settings/Distillation and cannot be changed.

Amount

Enter the amount, which the flask contains.

Distillate

Enter the amount, which you want to distill off.

Efficiency

Balancing requires a precise thermal efficiency to be entered. The efficiency is dependent on a combination of factors that affect the distillation process, so it must initially be estimated. The first time a distillation process is carried out should therefore be used for calibration. Following completion of the distillation process, determine the quantity of distillate produced and calculate the actual heat rejection rate according to the following formula:

$$\eta^p = \frac{\eta_{th} \cdot m_{(calc)}}{m_{(gem)}}$$

η^p	Actual efficiency
η_{th}	Estimated efficiency
$m_{(gem)}$	Measured quantity of distillate
$m_{(calc)}$	Specified (calculated) distillate quantity

For subsequent distillations, the actual calculated efficiency should be entered; this will ensure that subsequent automatic distillations under the same distillation and ambient conditions are performed with sufficient precision.

End of automatic distillation


Automatic distillation will stop when the following conditions are met:

- Volume controlled**
The specified quantity of distillate has been reached.
- 100% Distillation**
The measured cooling water temperature has fallen below a target value, indicating that the solvent is fully evaporated.

The distillation can also be stopped manually:


- Press the Rotating/ Pressing knob with the "Stop" display field selected.

"User-defined Distillation" menu

"User-defined Distillation" display screen

1. Operating mode display
2. Control field for distillation procedure
3. Vacuum display (actual value)
4. Vacuum display (target value)
5. Control field for „Start/ Stop“ distillation
6. Timer display
7. Vacuum-rotation-time graph
8. Display of flow rate, differential temperature, and condenser heat rejection rate
9. Rotary speed display (target value) and speed measurement unit
10. Error message in case of fault, otherwise status indicator


- ☞ Turn the Rotating/ Pressing knob to the right/left to select „Procedure“.
- ☞ Press the Rotating/ Pressing knob.
- ☞ The display will show the "Procedures" screen.
- ☞ Turn the Rotating/ Pressing knob to the right/left to select a procedure.
- ☞ Press the SET button.
- ☞ The procedure selected will be used for the subsequent distillation.

Note: The vacuum-rotation-time graph is displayed for each procedure. The parameters can be changed for each procedure. Please refer to the next section for instructions on how to do this. The sequence used for a manual distillation process can also be saved as a procedure. see "Manual Mode" section.

Changing a procedure

- ☞ Select the procedure to be changed and press the Rotating/ Pressing knob. The "Delete and create" display field will be activated.
- ☞ Press the Rotating/ Pressing knob with the "Delete and create" field selected to delete the procedure.

New vacuum and rotation values can now be entered to define a new step in procedure.

Vacuum

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Vacuum" from the on-screen options.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to the right/left to enter a value.
- ☞ Press the Rotating/ Pressing knob.

Duration

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Duration" from the on-screen options.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to the right/left to enter a value.
- ☞ Press the Rotating/ Pressing knob. This will define a new step in the procedure with the pair of values vacuum-duration and add it to the procedure. The diagram will be updated.

Rotation

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Rotation" from the on-screen options.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to the right/left to enter a value.
- ☞ Press the Rotating/ Pressing knob.

Duration

- ☞ Turn the Rotating/ Pressing knob to the right/left to select "Duration" from the on-screen options.
- ☞ Press the Rotating/ Pressing knob.
- ☞ Turn the Rotating/ Pressing knob to the right/left to enter a value.
- ☞ Press the Rotating/ Pressing knob. This will define a new step in the procedure with the pair of values rotation-duration and add it to the procedure. The diagram will be updated.

The vacuum-duration and rotation-duration steps can be repeated a maximum of twenty times.

SET key: Exit menu, store the newly created procedure and activate it as the current procedure.

ESC key: Exit menu and reset all changes.

Stopping user-defined distillation


The user-defined distillation will stop automatically once all steps in the currently active procedure have been completed. It is also possible to stop the process manually by pressing the Rotating/ Pressing knob with the "Stop" display field selected.

Stand-by operation

- ☞ Press the "Power" button.
- ☞ The device will enter stand-by.

- ☞ Press the "Power" button again.
- ⇒ The screen will display "Start" to indicate that the device is ready for use again.

Setting the lower end stop

Caution! Depending on the size of the piston, the setting angle of the rotation drive and the position of the heating bath and lift, the evaporator piston can be in contact with the bottom of the heating bath.

Caution - glass may break!

Limit the lower lift position with the variable end stop.

The lift can be moved to any desired position in manual mode using the "▲" and "▼" keys. There is no automatic cut-off if a collision occurs.

☞ Press the "▼" key until the lift reaches the desired position.

Note: The evaporator piston should be 2/3 immersed in the heating bath.

☞ To adjust the stop element (Q), press the centre button (R) on the front of the lift (Fig. 18).

☞ Move the stop (Q) to the desired position (Fig. 19)

☞ Press the "▲" key until the lift has reached the upper end stop.

Note: The path is limited from 0 – 6 cm.


Fig. 28

Fig. 29

Check that the limit stop is functioning properly:

☞ Lower the lift by holding down the "▼" key.

☞ The drive stops automatically when the preset desired lower end position has been reached.

☞ Move the drive back to the upper position.

To activate the safety lift following a long period of inactivity, use the motor to move the lift to the top or bottom position several times prior to starting distillation. (See the section entitled Safety Instructions - Safety Lift!)

Filling the evaporator piston

Prior to filling the evaporator piston, a vacuum controller is used to regulate the glass apparatus to the target pressure.

☞ Now fill the evaporator piston using the backfeed line, see fig. 30.

☞ Due to the vacuum present, the solvent is suctioned into the evaporator piston. This enables you to keep solvent loss due to suctioning to a minimum.

☞ You can also fill the evaporator piston manually prior to creating the vacuum. The evaporator piston should not be filled more than half its volume.

*Note: **Caution!** The maximum allowed load (evaporator piston and contents) is 3 kg.*

Fig. 30

Setting up the heating bath

Refer also to the operating instructions for the heating bath IKA® HB 10!

☞ Move the lift to the bottom position and check the position of the heating bath in relation to the evaporator piston. When using larger evaporator pistons (2 or 3 litres) or depending on the angle of the rotation drive, you can move the heating bath 50 mm to the right.

☞ Fill the heating bath with the tempering medium until the evaporator piston is surrounded by tempering medium to 2/3 of its volume.

☞ Switch on the rotation drive and slowly increase the speed

Note: Avoid creating waves.

☞ Switch on the heating bath using the main on/off switch.

Note: Avoid stress on the glass due to different evaporator piston and heating bath temperatures when lowering the evaporator piston into the heating bath!


Note: If non-original accessories are used that are not supplied by IKA®, then it is possible that the 50 mm travel range provided by the heating bath will not be sufficient. This applies in particular when using 3 L evaporator flasks with foam brake.

Use the IKA® RV 10.3000 mounting plate to extend the heating bath travel range by 150 mm.

Maintenance and cleaning

The device is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

The seal on the glass condenser should be checked at regular intervals and replaced if necessary.

 For cleaning disconnect the main plug.

To clean the device use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see www.ika.com, spare parts diagram and spare parts list.

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard. For this use the "**certificate of compliance**", download the form from the **IKA®** website www.ika.com.


If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Safety lift

The safety lift must be checked daily prior to operation!

After a long period of inactivity (approx. four weeks), use the motor to move the lift to the top or bottom position several times prior to starting distillation. (See the section entitled Safety Instructions - Safety Lift)! If the safety lift is not working, please contact the **IKA®** Service department.

Accessories

RV 10.1	NS 29/32 Vertical glassware (1)
RV 10.10	NS 29/32 Vertical glassware, coated (1)
RV 10.2	NS 29/32 Diagonal glassware (2)
RV 10.20	NS 29/32 Diagonal glassware, coated (2)

RV 10.700	NS 29/42 Vertical glassware (1)
RV 10.710	NS 29/42 Vertical glassware, coated (1)
RV 10.800	NS 29/42 Diagonal glassware (2)
RV 10.810	NS 29/42 Diagonal glassware, coated (2)

RV 10.900	NS 24/40 Vertical glassware (1)
RV 10.910	NS 24/40 Vertical glassware, coated (1)
RV 10.1000	NS 24/40 Diagonal glassware (2)
RV 10.1010	NS 24/40 Diagonal glassware, coated (2)

RV 10.3	Vertical-intensive condenser with manifold(1)
RV 10.30	Vertical-intensive condenser with manifold, coated (1)
RV 10.4	Dry ice condenser (2)
RV 10.40	Dry ice condenser, coated (2)
RV 10.5	Vertical-condenser with manifold and cut-off valve (no picture)
RV 10.50	Vertical-condenser with manifold and cut-off valve, coated (no picture)
RV 10.6	Vertical-intensive condenser with manifold and cut-off valve for reflux distillation (3)
RV 10.60	Vertical-intensive condenser with manifold and cut-off valve for reflux distillation (3)

HB 10.1	Protection shield (Heating bath HB 10) (1)
HB 10.2	Protection cover (Heating bath HB 10) (2)

RV 10.70	NS 29/32 Steam pipe
RV 10.71	NS 24/29 Steam pipe
RV 10.72	NS 29/42 Steam pipe
RV 10.73	NS 24/40 Steam pipe
RV 10.74	NS 29/32 Steam pipe, short, for reflux distillation
RV 10.75	NS 24/29 Steam pipe, short, for reflux distillation
RV 10.76	NS 29/42 Steam pipe, short, for reflux distillation
RV 10.77	NS 24/40 Steam pipe, short, for reflux distillation

RV 10.80	NS 29/32 Evaporator piston	50 ml
RV 10.81	NS 29/32 Evaporator piston	100 ml
RV 10.82	NS 29/32 Evaporator piston	250 ml
RV 10.83	NS 29/32 Evaporator piston	500 ml
RV 10.84	NS 29/32 Evaporator piston	1000 ml
RV 10.85	NS 29/32 Evaporator piston	2000 ml
RV 10.86	NS 29/32 Evaporator piston	3000 ml
RV 10.90	NS 24/32 Evaporator piston	50 ml
RV 10.91	NS 24/32 Evaporator piston	100 ml
RV 10.92	NS 24/32 Evaporator piston	250 ml
RV 10.93	NS 24/32 Evaporator piston	500 ml
RV 10.94	NS 24/32 Evaporator piston	1000 ml
RV 10.95	NS 24/32 Evaporator piston	2000 ml
RV 10.96	NS 24/32 Evaporator piston	3000 ml
RV 10.97	NS 24/40 Evaporator piston	1000 ml
RV 10.2001	NS 29/42 Evaporator piston	50 ml
RV 10.2002	NS 29/42 Evaporator piston	100 ml
RV 10.2003	NS 29/42 Evaporator piston	250 ml
RV 10.2004	NS 29/42 Evaporator piston	500 ml
RV 10.87	NS 29/42 Evaporator piston	1000 ml
RV 10.2005	NS 29/42 Evaporator piston	2000 ml
RV 10.2006	NS 29/42 Evaporator piston	3000 ml
RV 10.2007	NS 24/40 Evaporator piston	50 ml
RV 10.2008	NS 24/40 Evaporator piston	100 ml
RV 10.2009	NS 24/40 Evaporator piston	250 ml
RV 10.2010	NS 24/40 Evaporator piston	500 ml
RV 10.2011	NS 24/40 Evaporator piston	1000 ml
RV 10.2012	NS 24/40 Evaporator piston	2000 ml
RV 10.2013	NS 24/40 Evaporator piston	3000 ml

RV 10.100	KS 35/20 Receiving flask	100 ml
RV 10.101	KS 35/20 Receiving flask	250 ml
RV 10.102	KS 35/20 Receiving flask	500 ml
RV 10.103	KS 35/20 Receiving flask	1000 ml
RV 10.104	KS 35/20 Receiving flask	2000 ml
RV 10.105	KS 35/20 Receiving flask	3000 ml
RV 10.200	KS 35/20 Receiving flask, coated	100 ml
RV 10.201	KS 35/20 Receiving flask, coated	250 ml
RV 10.202	KS 35/20 Receiving flask, coated	500 ml
RV 10.203	KS 35/20 Receiving flask, coated	1000 ml
RV 10.204	KS 35/20 Receiving flask, coated	2000 ml
RV 10.205	KS 35/20 Receiving flask, coated	3000 ml

RV 10.300	NS 29/32 Powder piston	500 ml
RV 10.301	NS 29/32 Powder piston	1000 ml
RV 10.302	NS 29/32 Powder piston	2000 ml
RV 10.303	NS 24/29 Powder piston	500 ml
RV 10.304	NS 24/29 Powder piston	1000 ml
RV 10.305	NS 24/29 Powder piston	2000 ml
RV 10.2014	NS 29/32 Powder piston	500 ml
RV 10.2015	NS 29/32 Powder piston	1000 ml
RV 10.2016	NS 29/32 Powder piston	2000 ml
RV 10.217	NS 24/40 Powder piston	500 ml
RV 10.218	NS 24/40 Powder piston	1000 ml
RV 10.219	NS 24/40 Powder piston	2000 ml

RV 10.400	NS 29/32 Evaporator cylinder	500 ml
RV 10.401	NS 29/32 Evaporator cylinder	1500 ml

	RV 10.402 RV 10.403	NS 24/29 Evaporator cylinder 500 ml NS 24/29 Evaporator cylinder 1500 ml
	RV 10.2020 RV 10.2021	NS 29/42 Evaporator cylinder 500 ml NS 29/42 Evaporator cylinder 1500 ml
	RV 10.2022 RV 10.2023	NS 24/40 Evaporator cylinder 500 ml NS 24/40 Evaporator cylinder 1500 ml
	RV 10.500 RV 10.501 RV 10.2024 RV 10.2025	NS 29/32 Foam brake NS 24/29 Foam brake NS 29/42 Foam brake NS 24/40 Foam brake
	RV 10.600 RV 10.601 RV 10.602	NS 29/32 Distillation spider with 6 sleeves NS 29/32 Distillation spider with 12 sleeves NS 29/32 Distillation spider with 20 sleeves
	RV 10.603 RV 10.604 RV 10.605	NS 24/29 Distillation spider with 6 sleeves NS 24/29 Distillation spider with 12 sleeves NS 24/29 Distillation spider with 20 sleeves
	RV 10.2026 RV 10.2027 RV 10.2028	NS 29/42 Distillation spider with 6 sleeves NS 29/42 Distillation spider with 12 sleeves NS 29/42 Distillation spider with 20 sleeves
	RV 10.2029 RV 10.2030 RV 10.2031	NS 24/40 Distillation spider with 6 sleeves NS 24/40 Distillation spider with 12 sleeves NS 24/40 Distillation spider with 20 sleeves
	RV 10.610	Distilling sleeve 20 ml
	RV 10.606 RV 10.607	NS 29/32 Distillation spider with 5 flasks 50 ml NS 29/32 Distillation spider with 5 flasks 100 ml
	RV 10.608 RV 10.609	NS 24/29 Distillation spider with 5 flasks 50 ml NS 24/29 Distillation spider with 5 flasks 100 ml
	RV 10.2032 RV 10.2033	NS 29/42 Distillation spider with 5 flasks 50 ml NS 29/42 Distillation spider with 5 flasks 100 ml
	RV 10.2034 RV 10.2035	NS 24/40 Distillation spider with 5 flasks 50 ml NS 24/40 Distillation spider with 5 flasks 100 ml
	RV 10.3000	Extension plate
	RV 10.4001	Solenoid valve (central vacuum system) (1)
	RV 10.4002	Solenoid valve (laboratory vacuum system) (2)
	RV 10.4003 RV 10.5001	Pump controller, including solenoid valve (3) Water regulator valve (4), only for use on water supply line!
	RV 10.5002	Filter (5)
	RV 10.5003	Pressure control valve (6)
		
		
		
	RV 10.8001	Seal, PTFE compound
	PC 1.1	Cable Labworldsoft®

Error codes

Any malfunctions during operation will be identified by an error message on the display.

Once a serious error message has been displayed, the lift moves to the top end position and the device can no longer be operated. The lift can be operated again.

Proceed as follows in such cases:

- ☞ Switch off device using the main switch
- ☞ Carry out corrective measures
- ☞ Restart device

Error	Effect	Cause	Correction
Communication error bath	Distillation cancelled in "Auto-distillation" mode	IR interface obstructed Heating bath switched off	- Press ESC to confirm - Switch on heating bath - Check and clean IR interface
No rotation	Distillation stopped, no rotary drive	Rotary drive does not start or exhibits very large speed variation	- Switch device off and on at mains switch
No temperature difference	Distillation cancelled in "Auto-distillation" mode	No differential temperature increase of cooling water	- Press ESC to confirm - Check flow of cooling water - Check that the cooling water hoses are fitted correctly
Temperature out of range	Distillation stopped	Differential temperature outside measurement range (> 7 K)	- Press ESC to confirm - Check that the cooling water hoses are fitted correctly
No vacuum	Distillation stopped Vacuum cannot be controlled	Vacuum target value not reached	- Press ESC to confirm - Check that the glassware and hose connections are properly sealed - Switch on the pump - Check the solenoid valve connections
No venting	Vacuum cannot be vented	No venting of condenser	- Switch off the device - Vacuum sensor faulty - Manually vent the glassware using the stop valve
Water valve not closed	Water runs continuously into cooling circuit	Water valve not closed	- Press ESC to confirm - Check the connection
Flow rate out of range	Distillation cancelled in "Auto-distillation" mode	Flow rate of cooling water out of range	- Press ESC to confirm - Readjust flow rate - Switch off flow rate control (for dry ice coolers)
Condenser overload	Distillation stopped Distillate may be able to enter the vacuum pump	Condenser overload	- Press ESC to confirm - Change process parameters, e.g. - Increase pressure - Reduce bath temperature - Increase condenser flow rate within specified range
Distillation dry	Distillation stopped	The system detected that there was no more medium in the evaporating flask during the distillation.	- Press ESC to confirm - Alternatively, activate the "Drying" option in the "Settings" menu
Lift not in end position	Distillation stopped Lift motor switches off	Obstruction in path of lift	- Press ESC to confirm - Check the path of the lift. - Remove any obstructions.
Calibration error	Calibration reset to factory setting	Out of calibration range (max. +/-0.5 K)	- Press ESC to confirm - Check that the sensor is connected correctly - Carry out a new calibration procedure
Memory error	Calibration parameter for vacuum reset to factory setting	Error by checking contents of memory	Vacuum must be recalibrated
Vacuum leak	Distillation stopped	Leak in vacuum system	Press ESC to confirm Locate leak

Error	Effect	Cause	Correction
No boiling point found	Distillation stopped	The boiling point could not be determined	- High-boiling solvent: restart - Check heating bath
No analogue pump found	The distillation was not started	No analogue pump was found	- Confirm with ESC - Check pump and connections

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- contact the **IKA**® service department,
- send the device for repair, including a short description of the fault.

Technical Data

Operating voltage range	Vac	100 - 230 ±10%
Rated voltage	Vac	100 - 230
Frequency	Hz	50 / 60
Power input without heating bath	W	100
Power input operation "stand by"	W	3,3
Speed	rpm	20-280
Speed display		digital
Dimensions of visible display area (W x H)	mm	70 x 52
Display		TFT-Display
Multiple languages		yes
Right and left movement/ interval operation		yes
Smooth start		yes
Lift		automatic
Lifting speed	mm/s	50
Stroke	mm	140
Setting of lower end stop	mm	60, contactless
Head angle adjustable		0° - 45°
Timer	min	1-200
Interval	sec	1-100
Integrated vacuum controller		yes
Measurement range vacuum	mbar (hPa)	1050 - 1
Adjustment range vacuum	mbar (hPa)	1012 - 1
Measurement accuracy vacuum	mbar	± 2 (± 1 Digit) if calibration is done at a constant temperature
Adjustment accuracy vacuum	mbar	adjustable hysteresis
Measurement range of differential temperature	K	7
Offset balancing range	K	±0,5
Interface		USB, RS 232
Ramp programming		yes
Distillation processes programmable		yes
Remote control		with accessory labworldsoft®
Cooling surface	cm ²	1200
(Standard condenser RV 10.1,10.10,10.2,10.20)		
Minimum cooling water flow rate	l/h	30
Maximum cooling water flow rate	l/h	100
Cooling water pressure	bar	1
Cooling water temperature range	°C	18 - 22, constantly
Fuse		2xT1, 6A 250V 5x20
Perm. On- time	%	100
Perm. ambient temperature	°C	5-40
Perm. humidity	%	80
Protection acc. to DIN EN 60529		IP 20
Protection class		I
Overvoltage category		II
Contamination level		2
Weight (no glassware; no heating bath)	kg	18,6
Dimensions (W x D x H)	mm	500 x 440 x 430
Operation at a terrestrial altitude	m	max. 2000 above sea level

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GmbH & Co.KG