

reagents & CHEMICALS

agarose

Cleaver Scientific CleverGEL is an environmentally friendly agarose suitable for analysis of nucleic acids using standard electrophoretic procedures. Available as standard Low EEO, High Resolution PCR grade, Low melting point and Instant agarose tablets.



CleverGEL agarose is suitable for analysis of nucleic acids using standard electrophoretic procedures. It is manufactured by a process which excludes organic solvents harmful to marine life, making them far kinder to the environment than conventional agarose. A low EEO (electroendoosmotic) flow minimises diffusion so that even the smallest of nucleic acid fragments remains sharp and tightly resolved.

CleverGEL is available as standard Low EEO, High Resolution PCR grade, Low melting point and Instant agarose tablets.

KEY FEATURES

CleverGel Low EEO agarose:

- Ideal for routine analytical electrophoresis and blotting of DNA and RNA 0.1-10Kb in size
- Low EEO
- High gel strength

CleverGel High Resolution – PCR grade:

- High Resolution
- Low background for analysis of fragments 20-800bp

CleverGel Low Melting Point:

- Used for nucleic acid recovery
- typically resolves fragments 200bp to 25Kb

CleverGel Instant Agarose Tablets:

- Faster and simpler to prepare
- Low EEO
- Consistent gel percentage



Instant Agarose Tablets

TECHNICAL SPECIFICATIONS

	Low EEO	Low Melting Point	High Resolution	Instant Agarose
CAS	9012-36-6	39346-81-1	39346-81-1	9012-36-6
Gelling Point*	36°C±1.5°C	26-30°C	≤33°C	36°C±1.5°C
Melting Point*	88°C±1.5°C	≤65°C	≤70°C	88°C±1.5°C
Solubility	clear, colourless @ 1% [w/v] solution	clear, colourless @ 2% [w/v] solution		clear, colourless @ 1% [w/v] solution
Moisture	≤10%	≤10%	≤10%	≤10%
Gel Strength	>1200 g/cm ² (1% [w/v] gel)	>200 g/cm ² (1% [w/v] gel)	≥750 g/cm ² (1.5% [w/v] gel)	>1200 g/cm ² (1% [w/v] gel)
Nuclease & Protease Free	yes	yes	yes	yes

*For a 1.5% [w/v] gel

ORDERING INFORMATION

GENERAL PURPOSE

CSL-AG5	Agarose 5g, Low EEO
CSL-AG100	Agarose 100g, Low EEO
CSL-AG500	Agarose 500g, Low EEO
CSL-AG1000	Agarose 1000g, Low EEO (2x500g bottles)
CSL-AG2000	Agarose 2000g, Low EEO (4x500g)
CSL-AG5000	Agarose 5000g, Low EEO (10x500g)
CSL-AG10KG	Agarose 10Kg, Low EEO (20x500g)

LOW MELTING POINT

CSL-LMA5	Agarose 5g, LMP
CSL-LMA50	Agarose 50g, LMP
CSL-LMA100	Agarose 100g, LMP
CSL-HRA5	Agarose 5g, HR
CSL-HRA100	Agarose 100g, HR
CSL-HRA500	Agarose 500g, HR

AGAROSE TABLETS

CSL-AGTAB	Agarose 100g, Low EEO (200x 0.5g tablets, supplied as 20 blister packs of 10x 0.5g tablets)
-----------	---