




EASY, FAST, RELIABLE

ELDONCARD® 2521: Bed Side Test AB forward blood typing



Recommended for:

- Medical Market – Blood Banks, Hospitals, Clinics etc.
- Emergency & Rescue – NGO's
- Military, Police and Civil Service
- Donor Recruitment
- Rhesus tests for Pregnant Women
- Educational Market
- Health / Home Market
- Bed side test

ELDONCARD® 2521 carries 2 x 4 circles with a diameter of 20 mm. In the circles, dried blood typing reagents are found. The card is used for the detection of the blood type of one person (A, B, AB, or O and D (RhD positive or negative) in each of the 4 circles of one drop of water and one drop of blood from the testperson. Control field to detect false possible false positives. The procedure is described in detail in the instructions for use.

Available as:



- **ELDON BAG**
Contains:
 - 25 ELDONCARD 2521
 - 25 Blood Group Determinations



- **ELDON BOX 2521**
Contains:*
 - 100 ELDONCARDS 2521
 - 100 Blood Group Determinations

- **Costume made Kits for special use NGO/Doctors/ Military**

*Each EldonCard is in an aluminum envelope.

Labelling: CE & FDA approved and a large number of other countries

Storage: 5-37 Celsius / 41-99 Fahrenheit – **NO COOLCHAIN**

Expiration: 24 month from first day of production

Accessories:

In general, accessories are available in numbers matching the number of Blood Group Determinations.



- Safety Lancet



- Steel Lancet



- EldonStick



- Alcohol Swabs



- EldonFoil



- Pipette

General: To perform one Blood Group Determination 2511 you need 1 lancet, 4 EldonSticks, 1 EldonFoil, 1 pipette and 1 alcohol swab.

Contacts: Kasper Juul Hedegaard - kj@eldoncard.com - +4542138115